

METHOD #: 310.1 Approved for NPDES (Editorial Revision 1978)
TITLE: Alkalinity (Titrimetric, pH 4.5)
ANALYTE: Alkalinity
INSTRUMENTATION: Titration
STORET No. 00410

1.0 Scope and Application

- 1.1 This method is applicable to drinking, surface, and saline waters, domestic and industrial wastes.
- 1.2 The method is suitable for all concentration ranges of alkalinity; however, appropriate aliquots should be used to avoid a titration volume greater than 50 mL.
- 1.3 Automated titrimetric analysis is equivalent.

2.0 Summary of Method

- 2.1 An unaltered sample is titrated to an electrometrically determined end point of pH 4.5. The sample must not be filtered, diluted, concentrated, or altered in any way.

3.0 Comments

- 3.1 The sample should be refrigerated at 4°C and run as soon as practical. Do not open sample bottle before analysis.
- 3.2 Substances, such as salts of weak organic and inorganic acids present in large amounts, may cause interference in the electrometric pH measurements.
- 3.3 For samples having high concentrations of mineral acids, such as mine wastes and associated receiving waters, titrate to an electrometric endpoint of pH 3.9, using the procedure in: Annual Book of ASTM Standards, Part 31, "Water", p 115, D- 067, Method D, (1976).
- 3.4 Oil and grease, by coating the pH electrode, may also interfere, causing sluggish response.

4.0 Apparatus

- 4.1 pH meter or electrically operated titrator that uses a glass electrode and can be read to 0.05 pH units. Standardize and calibrate according to manufacturer's instructions. If automatic temperature compensation is not provided, make titration at 25 ±2 °C.
- 4.2 Use an appropriate sized vessel to keep the air space above the solution at a minimum. Use a rubber stopper fitted with holes for the glass electrode, reference electrode (or combination electrode) and buret.
- 4.3 Magnetic stirrer, pipets, flasks and other standard laboratory equipment.

4.4 Burets, Pyrex 50, 25 and 10 mL.

5.0 Reagents

- 5.1 Sodium carbonate solution, approximately 0.05 N: Place 2.5 ± 0.2 g (to nearest mg) Na_2CO_3 (dried at 250°C for 4 hours and cooled in desiccator) into a 1 liter volumetric flask and dilute to the mark.
- 5.2 Standard acid (sulfuric or hydrochloric), 0.1 N: Dilute 3.0 mL conc H_2SO_4 or 8.3 mL conc HCl to 1 liter with distilled water. Standardize versus 40.0 mL of 0.05 N Na_2CO_3 solution with about 60 mL distilled water by titrating potentiometrically to pH of about 5. Lift electrode and rinse into beaker. Boil solution gently for 3-5 minutes under a watch glass cover. Cool to room temperature. Rinse cover glass into beaker. Continue titration to the pH inflection point. Calculate normality using:

$$N = \frac{A \times B}{53.00 \times C}$$

where:

A = g Na_2CO_3 weighed into 1 liter

B = mL Na_2CO_3 solution

C = mL acid used to inflection point

- 5.3 Standard acid (sulfuric or hydrochloric), 0.02 N: Dilute 200.0 mL of 0.1000 N standard acid to 1 liter with distilled water. Standardize by potentiometric titration of 15.0 mL 0.05N Na_2CO_3 solution as above.

6.0 Procedure

6.1 Sample size

- 6.1.1 Use a sufficiently large volume of titrant (> 20 mL in a 50 mL buret) to obtain good precision while keeping volume low enough to permit sharp end point.
- 6.1.2 For < 1000 mg CaCO_3/L use 0.02 N titrant
- 6.1.3 For > 1000 mg CaCO_3/L use 0.1 N titrant
- 6.1.4 A preliminary titration is helpful.

6.2 Potentiometric titration

- 6.2.1 Place sample in flask by pipetting with pipet tip near bottom of flask
- 6.2.2 Measure pH of sample
- 6.2.3 Add standard acid (5.2 or 5.3), being careful to stir thoroughly but gently to allow needle to obtain equilibrium.
- 6.2.4 Titrate to pH 4.5. Record volume of titrant.

6.3 Potentiometric titration of low alkalinity

- 6.3.1 For alkalinity of < 20 mg/L titrate 100-200 mL as above (6.2) using a 10 mL microburet and 0.02 N acid solution (5.3).
- 6.3.2 Stop titration at pH in range of 4.3-4.7, record volume and exact pH. Very carefully add titrant to lower pH exactly 0.3 pH units and record volume.

7.0 Calculations

7.1 Potentiometric titration to pH 4.5

$$\text{Alkalinity, mg/L CaCO}_3 = \frac{A \times N \times 50,000}{\text{mL of sample}}$$

where:

A = mL standard acid

N = normality standard acid

7.2 Potentiometric titration of low alkalinity:

$$\text{Total Alkalinity, mg/L CaCO}_3 = \frac{(2B - C) \times N \times 50,000}{\text{mL of sample}}$$

where:

B = mL titrant to first recorded pH

C = total mL titrant to reach pH 0.3 units lower

N = normality of acid

8.0 Precision and Accuracy

8.1 Forty analysts in seventeen laboratories analyzed synthetic water samples containing increments of bicarbonate, with the following results:

Increment as Alkalinity mg/liter, CaCO ₃	Precision as Standard Deviation mg/liter, CaCO ₃	Accuracy as Bias, %	Bias, mg/L, CaCO ₃
8	1.27	+ 10.61	+0.85
9	1.14	+22.29	+2.0
113	5.28	- 8.19	-9.3
119	5.36	- 7.42	-8.8

(FWPCA Method Study 1, Mineral and Physical Analyses)

8.2 In a single laboratory (EMSL) using surface water samples at an average concentration of 122 mg CaCO₃/L, the standard deviation was ±3.

Bibliography

1. Standard Methods for the Examination of Water and Wastewater, 14th Edition, p 278, Method 403, (1975).
2. Annual Book of ASTM Standards, Part 31, "Water", p 113, D-1067, Method B, (1976).