



Orthophosphate, USEPA by Segmented Flow Analysis (SFA) or Flow Injection Analysis (FIA)


(Cartridge Part #A002596)

1.0 Scope and Application


- 1.1 This method is used for the determination of orthophosphate in drinking water, surface water, saline water, and domestic and industrial waste according to USEPA Method 365.1 (Reference 15.4).
- 1.2 The Method Detection Limit (MDL) of this method is 0.001 mg/L phosphorus (P). The applicable range of the method is 0.01–1.0 mg/L phosphorus. The range may be extended to analyze higher concentrations by sample dilution.

2.0 Summary of Method

- 2.1 Orthophosphate reacts with molybdenum(VI) and antimony(III) in an acidic solution to form an antimony-phosphomolybdate complex. This complex is subsequently reduced with ascorbic acid to form a blue color, and the absorbance is measured at 660 nm (Reference 15.4).
- 2.2 The quality of the analysis is assured through reproducible calibration and testing of the Segmented Flow Analysis (SFA) or Flow Injection Analysis (FIA) system.
- 2.3 A general flow diagram of the SFA system is shown below (see Section 17.0 for a detailed flow diagram).



- 2.4 A general flow diagram of the FIA system is shown below (see Section 17.0 for a detailed flow diagram).



3.0 Definitions

Definitions for terms used in this method are provided in Section 16.0, "Glossary of Definitions and Purposes."

4.0 Interferences

- 4.1 Filter turbid samples prior to analysis.
- 4.2 The presence of less than 50 mg/L of iron(III), less than 10 mg/L of copper, or less than 10 mg/L of silicate does not interfere with this assay.
- 4.3 Samples with background absorbance at the analytical wavelength may interfere.

5.0 Safety

- 5.1 The toxicity or carcinogenicity of each compound or reagent used in this method has not been fully established. Each chemical should be treated as a potential health hazard. Exposure to these chemicals should be reduced to the lowest possible level.
- 5.2 For reference purposes, a file of Material Safety Data Sheets (MSDS) for each chemical used in this method should be available to all personnel involved in this chemical analysis. The preparation of a formal safety plan is also advisable.
- 5.3 The following chemicals used in this method may be highly toxic or hazardous and should be handled with extreme caution at all times. Consult the appropriate MSDS before handling.

5.3.1 Ammonium Molybdate Tetrahydrate, $(\text{NH}_4)_6\text{Mo}_7\text{O}_{24} \cdot 4\text{H}_2\text{O}$ (FW 1,235.95)

5.3.2 Antimony Potassium Tartrate Hemihydrate, $\text{K}(\text{SbO})\text{C}_4\text{H}_4\text{O}_6 \cdot \frac{1}{2}\text{H}_2\text{O}$ (FW 324.92)

- 5.3.3 Hydrochloric Acid, concentrated, HCl (FW 36.46)
- 5.3.4 Potassium Phosphate Monobasic, KH_2PO_4 (FW 136.09)
- 5.3.5 Sodium Hydroxide, NaOH (FW 40.00)
- 5.3.6 Sulfuric Acid, concentrated, H_2SO_4 (FW 98.08)
- 5.4 Unknown samples may be potentially hazardous and should be handled with extreme caution at all times.
- 5.5 Proper personal protective equipment (PPE) should be used when handling or working in the presence of chemicals.
- 5.6 This method does not address all safety issues associated with its use. The laboratory is responsible for maintaining a safe work environment and a current awareness file of OSHA regulations regarding the safe handling of the chemicals specified in this method.

6.0 Apparatus, Equipment, and Supplies

- 6.1 Segmented Flow Analysis (SFA) System (OI Analytical Flow Solution[®] IV) consisting of the following:
 - 6.1.1 Model 502 Multichannel Peristaltic Pump
 - 6.1.2 Random Access (RA) Autosampler
 - 6.1.3 Expanded Range (ER) Photometric Detector with 5-mm path length flowcell and 660-nm optical filter
 - 6.1.4 Data Acquisition System (PC or Notebook PC) with WinFLOW[™] software
 - 6.1.5 Orthophosphate/Total Phosphorus, USEPA Cartridge (Part #A002596)
 - 6.1.6 For FIA, Flow Solution IV must be equipped with the FIA option.
- 6.2 Sampling equipment—Sample bottle, amber glass, with polytetrafluoroethylene (PTFE)-lined cap. Clean by washing with detergent and water, rinsing with two aliquots of reagent water, and drying by baking at 110°–150°C for a minimum of one hour.
- 6.3 Standard laboratory equipment including volumetric flasks, pipettes, syringes, etc. should all be cleaned, rinsed, and dried per bottle cleaning procedure in Section 6.2.

7.0 Reagents and Calibrants

7.1 Raw Materials

- 7.1.1 Ammonium Molybdate Tetrahydrate, $(\text{NH}_4)_6\text{Mo}_7\text{O}_{24}\cdot 4\text{H}_2\text{O}$ (FW 1,235.95)
- 7.1.2 Antimony Potassium Tartrate Hemihydrate, $\text{K}(\text{SbO})\text{C}_4\text{H}_4\text{O}_6\cdot \frac{1}{2}\text{H}_2\text{O}$ (FW 324.92)
- 7.1.3 Ascorbic Acid, $\text{C}_6\text{H}_8\text{O}_6$ (FW 176.12)
- 7.1.4 Deionized Water (ASTM Type I or II)
- 7.1.5 DOWFAX® 2A1 (Part #A000080)
- 7.1.6 Hydrochloric Acid, concentrated, HCl (FW 36.46)
- 7.1.7 Potassium Phosphate Monobasic, KH_2PO_4 (FW 136.09)
- 7.1.8 Sodium Hydroxide, NaOH (FW 40.00)
- 7.1.9 Sulfuric Acid, concentrated, H_2SO_4 (FW 98.08)

7.2 Reagent Preparation

Note: For best results, filter and degas all reagents prior to use.

7.2.1 Reagent Water

7.2.1.1 Degassed and deionized reagent water can be prepared in one of the following manners:

7.2.1.1.1 Place distilled/deionized water under a strong vacuum for 15–20 minutes. Magnetic stirring or sonification will aid in the degassing process.

7.2.1.1.2 Purge distilled/deionized water with a stream of nitrogen gas (or other inert gas) through a glass frit for approximately 5 minutes.

7.2.1.1.3 Boil distilled/deionized water in an Erlenmeyer flask for 15–20 minutes. Remove the flask from the heat source, cover it with an inverted beaker, and allow it to cool to room temperature.

7.2.1.2 After preparing the degassed reagent water, store the reagent water in a tightly sealed container to protect it from reabsorption of atmospheric gases. For best results, store degassed reagent water under a slight vacuum when not in use.

7.2.2 Start-up Solution/Diluent (1 L)

7.2.2.1 Add 2 mL of DOWFAX 2A1 to approximately 800 mL of reagent water (Section 7.2.1) in a 1-L volumetric flask.

7.2.2.2 Dilute to 1,000 mL with reagent water and mix gently.

7.2.3 Stock Sulfuric Acid, 5 N (1 L)

7.2.3.1 While stirring, carefully add 140 mL of concentrated sulfuric acid to approximately 700 mL of reagent water in a 1-L volumetric flask.

7.2.3.2 Allow to cool to room temperature. Dilute to 1,000 mL with reagent water and mix well.

Warning: Mixing sulfuric acid with water releases a great amount of heat. Take appropriate precautions.

7.2.4 Stock Antimony Potassium Tartrate Solution (1 L)

7.2.4.1 Dissolve 3 g of antimony potassium tartrate hemihydrate in approximately 800 mL of reagent water in a 1-L volumetric flask.

7.2.4.2 Dilute to 1,000 mL with reagent water and mix well.

Note: Store in an amber bottle at 4°C. This reagent is stable for 4–6 weeks if stored properly.

7.2.5 Stock Ammonium Molybdate Solution (1 L)

7.2.5.1 Dissolve 40 g of ammonium molybdate tetrahydrate in approximately 800 mL of reagent water in a 1-L volumetric flask.

7.2.5.2 Dilute to 1,000 mL with reagent water and mix well.

Note: Store in a polyethylene bottle at 4°C. This reagent is stable for 4–6 weeks if stored properly.

7.2.6 Stock Ascorbic Acid Solution (1 L)

7.2.6.1 Dissolve 18 g of ascorbic acid in approximately 800 mL of reagent water in a 1-L volumetric flask.

7.2.6.2 Dilute to 1,000 mL with reagent water and mix well.

Note: Store in an amber bottle at 4°C. This reagent is stable for 1 week if stored properly.

7.2.7 Color Reagent (200 mL)

7.2.7.1 Add the following solutions to a 200-mL volumetric flask in the order stated. Mix after each addition to prevent the solution from darkening.

7.2.7.1.1 Stock Sulfuric Acid (Section 7.2.3)—100 mL

7.2.7.1.2 Stock Antimony Potassium Tartrate Solution (Section 7.2.4)—10 mL

7.2.7.1.3 Stock Ammonium Molybdate Solution (Section 7.2.5)—30 mL

7.2.7.1.4 Stock Ascorbic Acid Solution (Section 7.2.6)—60 mL

7.2.7.1 Add 0.5 mL of DOWFAX 2A1 and mix gently.

Note: Prepare this reagent fresh daily.

7.2.8 0.1 N Hydrochloric Acid (500 mL)

7.2.8.1 While stirring, carefully add 4.15 mL of concentrated hydrochloric acid to approximately 400 mL of reagent water in a 500-mL volumetric flask.

7.2.8.2 Dilute to 500 mL with reagent water and mix well.

Warning: Mixing hydrochloric acid with water releases a great amount of heat. Take appropriate precautions.

7.2.9 1 N Sodium Hydroxide (500 mL)

7.2.9.1 While stirring, carefully add 20 g of sodium hydroxide to approximately 400 mL of reagent water in a 500-mL volumetric flask.

7.2.9.2 Cool the solution to room temperature. Dilute to 500 mL with reagent water and mix well.

Warning: Mixing sodium hydroxide with water releases a great amount of heat. Take appropriate precautions.

7.2.10 Carrier—Reagent Water

7.3 Calibrant Preparation

7.3.1 Stock Calibrant 1,000 mg/L Phosphorus (1 L)

7.3.1.1 Dissolve 4.393 g of potassium phosphate monobasic (dried at 100°C for 2 hours) in approximately 800 mL of reagent water in a 1-L volumetric flask.

7.3.1.2 Dilute to 1,000 mL with reagent water and mix well.

Note: Store in at 4°C. This reagent is stable for 4–6 weeks.

7.3.2 Intermediate Calibrant 100 mg/L Phosphorus (100 mL)

7.3.2.1 Use a volumetric pipet to add 10 mL of stock calibrant (Section 7.3.1) to approximately 80 mL of reagent water in a 100-mL volumetric flask.

7.3.2.2 Dilute to 100 mL with reagent water and mix well.

Note: Prepare the intermediate calibrant daily.

7.3.3 Working Calibrants (100 mL)

7.3.3.1 Add the designated volumes of stock calibrant (see Equation 1) to the required number of 100-mL volumetric flasks that each contain approximately 80 mL of reagent water.

7.3.3.2 Dilute each solution to the mark with reagent water and mix well.

Note: Prepare the working calibrants daily.

EQUATION 1

$$C_1 V_1 = C_2 V_2$$

Where:

C_1 = Concentration (in mg/L) of stock solution (or calibrant)

V_1 = Volume (in L) of stock solution (or calibrant) to be used

C_2 = Desired concentration (in mg/L) of working calibrant to be prepared

V_2 = Final volume (in L) of working calibrant to be prepared

By solving this equation for the volume of stock solution to be used (V_1), the following equation is obtained:

$$V_1 = \frac{C_2 V_2}{C_1}$$

Since the desired concentration (C_2), the final volume (V_2), and the concentration of the stock solution (C_1) are all known for any given calibrant concentration in a defined volume, the volume of stock solution to be used (V_1) is easily calculated.

7.3.3.3 Calibrants covering the entire range of this analysis can be prepared from the following table.

Final Concentration (mg/L)	Vol. of Inter. Cal. (μL)	Conc. of Inter. Cal. (mg/L)	Final Volume (mL)
0.01	10	100	100
0.20	200	100	100
0.40	400	100	100
0.60	600	100	100
0.80	800	100	100
1.0	1,000	100	100
2.0	2,000	100	100

8.0 Sample Collection, Preservation, and Storage

- 8.1 Samples should be collected in plastic or glass bottles that have been thoroughly cleaned and rinsed with reagent water (Section 7.2.1).
- 8.2 The volume of sample collected should be sufficient to ensure that a representative sample is obtained, replicate analysis is possible, and waste disposal is minimized.
- 8.3 Determine orthophosphate in unpreserved samples immediately upon collection.
- 8.4 Preserve and store samples at 4°C. Sample analysis should be performed as soon as possible to eliminate loss of analyte.
- 8.5 Holding time for preserved samples is 48 hours from the time of collection (Reference 15.5).

9.0 Quality Control

- 9.1 Each laboratory that uses this method is required to operate a formal quality assurance program (Reference 15.2). The minimum requirements of this program consist of an initial demonstration of laboratory capability and the periodic analysis of Laboratory Control Samples (LCSs) and Matrix Spike/Matrix Spike Duplicates (MS/MSDs) as a continuing check on performance. Laboratory performance is compared to established performance criteria to determine if the results of the analyses meet the performance characteristics of the method.
 - 9.1.1 The analyst shall make an initial demonstration of the ability to generate acceptable precision and accuracy with this method. This ability is established as described in Section 9.2.
 - 9.1.2 In recognition of advances that are occurring in analytical technology and to allow the analyst to overcome sample matrix interferences, the analyst is permitted certain options to improve performance or lower the costs of measurements. Alternate determinative techniques, such as the substitution of spectroscopic or other techniques, and changes that degrade method performance are not allowed. If an analytical technique other than the techniques specified in this method is used, that technique must have a specificity equal to or better than the specificity of the techniques in this method for the analyte(s) of interest.
 - 9.1.2.1 Each time a modification is made to this method, the analyst is required to repeat the procedure in Section 9.2. If the detection limit of the method will be affected by the change, the laboratory is required to demonstrate that the MDL is lower than one-third the regulatory compliance level or as low as or lower than that listed in Section 1.2. If calibration will be affected by the change, the analyst must recalibrate the instrument per Section 10.4.
 - 9.1.2.2 The laboratory is required to maintain records of modifications made to this method. These records include the information in this subsection, at a minimum.

- 9.1.2.2.1 The names, titles, addresses, and telephone numbers of the analyst(s) who performed the analyses and modification, and of the quality control officer who witnessed and will verify the analyses and modification.
- 9.1.2.2.2 A narrative stating the reason(s) for the modification.
- 9.1.2.2.3 Results from all quality control (QC) tests comparing the modified method to this method including:
- a) calibration (Section 10.4)
 - b) calibration verification (Section 9.5)
 - c) initial precision and recovery (Section 9.2.2)
 - d) analysis of blanks (Section 9.4)
 - e) ongoing precision and recovery (Section 9.6)
 - f) matrix spike and matrix spike duplicate (Section 9.3)
- 9.1.2.2.4 Data that will allow an independent reviewer to validate each determination by tracing the instrument output (peak height, area, or other signal) to the final result. These data are to include:
- a) sample numbers and other identifiers
 - b) analysis dates and times
 - c) analysis sequence/run chronology
 - d) sample weight or volume
 - e) sample volume prior to each cleanup step, if applicable
 - f) sample volume after each cleanup step, if applicable
 - g) final sample volume prior to injection
 - h) injection volume
 - i) dilution data, differentiating between dilution of a sample or modified sample
 - j) instrument and operating conditions
 - k) other operating conditions

- l) detector
 - m) printer tapes, disks, and other recording of raw data
 - n) quantitation reports, data system outputs, and other data necessary to link raw data to the results reported
- 9.1.3 Analyses of MS/MSD samples are required to demonstrate method accuracy and precision and to monitor matrix interferences (interferences caused by the sample matrix). The procedure and QC criteria for spiking are described in Section 9.3.
- 9.1.4 Analyses of laboratory reagent blanks (LRBs) are required to demonstrate freedom from contamination and that the compounds of interest and interfering compounds have not been carried over from a previous analysis. The procedures and criteria for analysis of an LRB are described in Section 9.4.
- 9.1.5 The laboratory shall, on an ongoing basis, demonstrate through the analysis of the LCS that the analytical system is in control. This procedure is described in Section 9.6.
- 9.1.6 The laboratory should maintain records to define the quality of data that is generated. Development of accuracy statements is described in Sections 9.3.8 and 9.6.3.
- 9.1.7 Accompanying QC for the determination of orthophosphate is required per analytical batch. An analytical batch is a set of samples analyzed at the same time to a maximum of 10 samples. Each analytical batch of 10 or fewer samples must be accompanied by a laboratory reagent blank (LRB, Section 9.4), a laboratory control sample (LCS, Section 9.6), and a matrix spike and matrix spike duplicate (MS/MSD, Section 9.3), resulting in a minimum of five analyses (1 sample, 1 LRB, 1 LCS, 1 MS, and 1 MSD) and a maximum of 14 analyses (10 samples, 1 LRB, 1 LCS, 1 MS, and 1 MSD) in the batch. If more than 10 samples are analyzed at one time, the samples must be separated into analytical batches of 10 or fewer samples.
- 9.2 Initial Demonstration of Laboratory Capability
- 9.2.1 Method Detection Limit (MDL)—To establish the ability to detect orthophosphate at low levels, the analyst shall determine the MDL per the procedure in 40 CFR 136, Appendix B (Reference 15.1) using the apparatus, reagents, and standards that will be used in the practice of this method. An MDL less than or equal to the MDL listed in Section 1.2 must be achieved prior to practice of this method.
- 9.2.2 Initial Precision and Recovery (IPR)—To establish the ability to generate acceptable precision and accuracy, the analyst shall perform the following operations:
- 9.2.2.1 Analyze four samples of the LCS (Section 9.6) according to the procedure beginning in Section 10.0.
 - 9.2.2.2 Using the results of the set of the four analyses, compute the average percent recovery (\bar{x}) and the standard deviation of the percent recovery (s) for orthophosphate. Use Equation 2 for the calculation of the standard deviation of the percent recovery (s).

EQUATION 2

$$s = \sqrt{\frac{\sum x^2 - \frac{(\sum x)^2}{n}}{n-1}}$$

Where:

s = Standard deviation

n = Number of samples

x = Percent recovery in each sample

9.2.2.3 Compare s and x with the precision and percent recovery acceptance criteria specified in Section 13.0. If the value of s exceeds the precision limit or the value of x falls outside the range for recovery, system performance is unacceptable and the problem must be found and corrected before the analysis may continue.

9.3 Matrix Spike/Matrix Spike Duplicate (MS/MSD)—The laboratory shall spike, in duplicate, a minimum of 10% of all samples (one sample in duplicate in each batch of 10 samples) from a given sampling site.

9.3.1 The concentration of the spike in the sample shall be determined as follows:

9.3.1.1 If, as in compliance monitoring, the concentration of orthophosphate in the sample is being checked against a regulatory concentration limit, the spiking level shall be at that limit or at one to five times higher than the background concentration of the sample (determined in Section 9.3.2), whichever concentration is higher.

9.3.1.2 If the concentration of orthophosphate in a sample is not being checked against a limit, the spike shall be at the concentration of the LCS or at least four times greater than the MDL.

9.3.2 Analyze one sample aliquot out of each set of 10 samples from each site or discharge according to the procedure beginning in Section 10.0 to determine the background concentration of orthophosphate.

9.3.2.1 If necessary, prepare a stock solution appropriate to produce a concentration level in the sample at the regulatory compliance limit or at one to five times the background concentration of orthophosphate (Section 9.3.1).

9.3.2.2 Spike two additional sample aliquots with the spiking solution (Section 9.3.2.1) and analyze these aliquots to determine the concentration after spiking.

9.3.3 Calculate the percent recovery of orthophosphate in each aliquot using Equation 3.

EQUATION 3

$$P = \frac{A - B}{T} \times 100$$

Where:

P = Percent recovery

A = Measured concentration of orthophosphate after spiking (Section 9.3.2.2)

B = Measured background concentration of orthophosphate (Section 9.3.2)

T = True concentration of the spike

9.3.4 Compare the recovery to the QC acceptance criteria in Section 13.0. If percent recovery is outside of the acceptance criteria, and the recovery of the LCS in the ongoing precision and recovery test (Section 9.6) for the analytical batch is within the acceptance criteria, an interference is present. In this case, the result may not be reported for regulatory compliance purposes.

9.3.5 If the results of both the MS/MSD and the LCS test fail the acceptance criteria, the analytical system is judged to be out of control. In this case, the problem shall be identified and corrected, and the analytical batch must be reanalyzed.

9.3.6 Compute the relative percent difference (RPD) between the two spiked sample results (Section 9.3.2.2, not between the two percent recoveries) using Equation 4.

EQUATION 4

$$RPD = \left[\frac{|D_1 - D_2|}{(D_1 + D_2) / 2} \right] \times 100$$

Where:

RPD = Relative percent difference

*D*₁ = Concentration of orthophosphate in the spiked sample

*D*₂ = Concentration of orthophosphate in the spiked duplicate sample

- 9.3.7 If the RPD is greater than 10%, the analytical system is judged to be out of control, and the problem must be immediately identified and corrected. The analytical batch must be reanalyzed.
- 9.3.8 As part of the QC program for the laboratory, method precision and accuracy for samples should be assessed and records should be maintained. After the analysis of five spiked samples in which the recovery passes the test in Section 9.3.4, compute the average percent recovery (P_a) and the standard deviation of the percent recovery (s_p). Express the accuracy assessment as a percent recovery interval from $P_a - 2s_p$ to $P_a + 2s_p$. For example, if $P_a = 90\%$ and $s_p = 10\%$ for five analyses, the accuracy interval is expressed as 70–110%. Update the accuracy assessment on a regular basis (e.g., after each 5–10 new accuracy measurements).
- 9.4 Laboratory Reagent Blanks (LRB)—Laboratory reagent blanks are analyzed to demonstrate freedom from contamination.
- 9.4.1 Analyze an LRB initially (i.e., with the tests in Section 9.2) and with each analytical batch. The LRB must be subjected to the exact same procedural steps as a sample.
- 9.4.2 If orthophosphate is detected in the LRB at a concentration greater than the ML, analysis of samples is halted until the source of contamination is eliminated and consequent analysis of another LRB shows no evidence of contamination.
- 9.5 Calibration Verification—Verify calibration of the analytical equipment before and after each analytical batch of 14 or fewer measurements. (The 14 measurements will normally be 10 samples, 1 LRB, 1 LCS, 1 MS, and 1 MSD). This can be accomplished by analyzing the midrange calibration standard and verifying that it is within the QC acceptance criteria for recovery in Section 13.0. (The concentration of the calibration verification depends on the calibration range being used.) Failure to attain recoveries within the acceptance criteria requires recalibration of the analytical system (Section 10.4).
- 9.6 Laboratory Control Sample (LCS)—To demonstrate that the analytical system is in control and acceptable precision and accuracy is being maintained with each analytical batch, the analyst shall perform the following operations:
- 9.6.1 Analyze an LCS with each analytical batch according to the procedure in Section 10.0.
- 9.6.2 If the precision and recovery for the LCS are within the acceptance criteria specified in Section 13.0, analysis of the batch may continue. If, however, the concentration is not within this range, the analytical process is not in control. In this event, correct the problem, repeat the LCS test, and reanalyze the batch.
- 9.6.3 The laboratory should add results that pass the specification in Section 9.6.2 to IPR and previous LCS data and update QC charts to form a graphic representation of continued laboratory performance. The laboratory should also develop a statement of laboratory data quality for orthophosphate by calculating the average percent recovery (R) and the standard deviation of the percent recovery (s_r). Express the accuracy as a recovery interval from $R - 2s_r$ to $R + 2s_r$. For example, if $R = 95\%$ and $s_r = 5\%$, the accuracy is 85–105%.

- 9.7 Reference Sample—To demonstrate that the analytical system is in control, the laboratory may wish to periodically test an external reference sample, such as a Standard Reference Material (SRM) available from the National Institute of Standards and Technology (NIST). Corrective action should be taken if the measured concentration significantly differs from the stated concentration.

10.0 Configuration and Start-up

10.1 Instrument Configuration

- 10.1.1 Configure the OI Analytical Flow Solution IV Analyzer according to the Operator's Manual and verify that each module is properly powered on.
- 10.1.2 Verify that the Orthophosphate/Total Phosphorus, USEPA Cartridge (Part #A002596) is configured as illustrated in the flow diagram shown in Section 17.0.
- 10.1.3 Connect the appropriate pump tubes to the cartridge and to their appropriate reagent containers according to the flow diagram.

10.2 Instrument Stabilization

- 10.2.1 Connect the reagent pump tubes to a reagent bottle containing the start-up solution (Section 7.2.2). Start the pump at 40% speed, allowing the start-up solution to flow through the entire system.
- 10.2.2 Verify that the flowcell of each detector is purged of all bubbles and the flow is stable and free from surging before proceeding.

10.3 Baseline Verification

- 10.3.1 Create and save a Method in WinFLOW. Refer to the WinFLOW Operator's Manual (Reference 15.7) for help on creating a Method.
- 10.3.2 Create and save a Sample Table in WinFLOW that will be used to generate a calibration curve using at least three calibrants that cover the full range of expected concentrations in the samples to be analyzed. This Sample Table should also be used to analyze all necessary QC samples as well as the analytical batch of samples to be analyzed. For help on creating a Sample Table, refer to the WinFLOW Operator's Manual (Reference 15.7).
- 10.3.3 Select **Collect Data** in the WinFLOW main window, enter the user's identification, select the appropriate Method and Sample Table, and begin to collect baseline data. Very sharp fluctuations in the baseline and/or consistent drifting are typically signs of bubbles in the flowcell. The flowcell must be free of bubbles prior to beginning analysis.

10.4 Calibration and Standardization

- 10.4.1 Prepare a series of at least three working calibrants using the stock solutions (Section 7.3) according to Equation 1, covering the desired analysis range.

- 10.4.2 Place the calibrants in the autosampler in order of decreasing concentration. Each calibrant should be analyzed according to the analytical procedures in Section 11.0. A calibration curve will be calculated by the WinFLOW software.
- 10.4.3 Acceptance or control limits for the calibration results should be established using the difference between the measured value of each calibrant and the corresponding “true” concentration.
- 10.4.4 Each calibration curve should be verified by analysis of a Laboratory Control Sample (LCS, Section 9.5). Using WinFLOW software, calibration, verification, and sample analysis may be performed in one continuous analysis.

11.0 Procedure

11.1 Analysis

- 11.1.1 Set the temperature of the heater to 37°C. Begin pump flow with the start-up solution (Section 7.2.2). Once the heater unit has reached 37°C, verify a stable baseline (Section 10.3).
- 11.1.2 After the baseline has been verified, place all reagents on-line and allow to pump at least 10–15 minutes and verify there are no bubbles in the flowcell. Obtain a stable baseline at 660 nm and autozero the baseline before beginning the analysis.
- 11.1.3 Load the sampler tray with calibrants, blanks, samples, and QC samples.

Note: The matrix of the working standards, blanks, and QC samples should match that of the samples being analyzed.
- 11.1.4 Using the Method and Sample Table created for the analytical batch to be analyzed and with the baseline verified to be stable, begin the analysis by selecting the “Fast Forward” button on the left side of the Data Analysis window in WinFLOW. This will initiate the sequential analysis of samples as defined in the Sample Table.
- 11.1.5 When analysis is complete, pump start-up solution through the system for at least 10–15 minutes. Stop the pump, release the tension on all pump tubes, and power off the system.

11.2 Operating Notes

- 11.2.1 For FIA operation, configure the cartridge as illustrated in Section 17.0, Figure 2.
- 11.2.2 The stock ascorbic acid solution (Section 7.2.6) should be a clear liquid, colorless to light amber in color. A darker color indicates poor reagent quality or an expired solution. Either remake the solution or change to a new source of ascorbic acid.
- 11.2.3 When starting up the system, run start-up solution (Section 7.2.2) for 10–15 minutes. Add the color reagent (Section 7.2.7) and wait for a stable baseline. The pH from the flowcell waste line should be approximately 1 when checked with pH paper.

11.2.4 If the baseline is noisy:

11.2.4.1 Check the color of the stock and working ascorbic acid solutions.

11.2.4.2 Filter all reagents prior to use through a 0.45- μ m membrane filter.

11.2.4.3 Acid wash sample cups, reagent bottles, and the wash reservoir with 0.1 N hydrochloric acid (Section 7.2.8) before use.

11.2.4.4 Detach the flowcell from the debubbler. Flush the flowcell forcefully with start-up solution to dislodge bubbles.

11.2.5 If abnormal or short peaks are observed:

11.2.5.1 Check that the solutions are flowing freely and consistently through the pump tubes.

11.2.5.2 Check that the temperature setting is correct.

11.2.6 To prevent carryover, wash the entire system with 1 N sodium hydroxide (Section 7.2.9) for 5 minutes. For consistent results, repeat this procedure every day before use and as needed. For persistent carryover problems, increase the sampler wash time.

11.2.7 Clear persistent air bubbles in the flowcell by creating back pressure. Pinch off the waste line from the flow cell for several seconds, then quickly release the line. Repeat as necessary.

11.2.8 If air spikes occur frequently:

11.2.8.1 Clean the system with 1 N sodium hydroxide (Section 7.2.9).

11.2.8.2 Change the debubbler pull-off pump tube to the next larger size.

11.2.8.3 Detach the flowcell from the debubbler and clean the flowcell with start-up solution (Section 7.2.2).

12.0 Data Analysis and Calculations

12.1 The calibration curve allows for accurate quantitation of the concentration in each sample.

12.2 WinFLOW software reports the concentration of each sample relative to the calibration curve.

13.0 Method Performance

Range:	0.01–1.0 mg/L
Throughput:	80 samples/hour
Precision:	
0.2 mg/L	<2% RSD
0.8 mg/L	<1% RSD
Method Detection Limit (MDL):	0.001 mg/L

14.0 Pollution Prevention and Waste Management

14.1 It is the laboratory's responsibility to comply with all federal, state, and local regulations governing waste management, particularly the hazardous waste identification rules and land-disposal restrictions. In addition, it is the laboratory's responsibility to protect air, water, and land resources by minimizing and controlling all releases from fume hoods and bench operations. Also, compliance is required with any sewage discharge permits and regulations.

14.2 For further information on waste management, consult Section 13.6 of *Less is Better: Laboratory Chemical Management for Waste Reduction* (Reference 15.3).

15.0 References

- 15.1 *Code of Federal Regulations*, Part 136, Title 40, Appendix B, 1994.
- 15.2 *Handbook for Analytical Quality Control in Water and Wastewater Laboratories*; EPA-600/4-79-019; U.S. Environmental Protection Agency, Office of Research and Development, Environmental Monitoring and Support Laboratory: Cincinnati, OH, 1979.
- 15.3 *Less is Better: Laboratory Chemical Management for Waste Reduction*. Available from the American Chemical Society, Department of Government Regulations and Science Policy, 1155 16th Street, NW, Washington, DC, 20036.
- 15.4 Phosphorus, All Forms. *Methods for Chemical Analysis of Water and Wastewater*; EPA-600/4-79-020; U.S. Environmental Protection Agency, Office of Research and Development, Environmental Monitoring and Support Laboratory: Cincinnati, OH, 1984; Method 365.1.
- 15.5 Sample Preservation. *Methods for Chemical Analysis of Water and Wastes*; EPA-600/4-79-020; U.S. Environmental Protection Agency, Office of Research and Development, Environmental Monitoring and Support Laboratory: Cincinnati, OH, 1984; xvii.
- 15.6 *Standard Methods for the Examination of Water and Wastewater*, 20th ed.; American Public Health Association: Washington, D.C., 1998.
- 15.7 WinFLOW Software and Operator's Manual (Part #A002877). Available from OI Analytical, P.O. Box 9010, College Station, TX, 77842-9010.

16.0 Glossary of Definitions and Purposes

The definitions and purposes are specific to this method but have been conformed to common usage as much as possible.

16.1 Units of weights and measures and their abbreviations

16.1.1 Symbols

°C	degrees Celsius
%	percent
±	plus or minus
≥	greater than or equal to
≤	less than or equal to

16.1.2 Alphabetical characters

g	gram
L	liter
mg	milligram
mg/L	milligram per liter
µg	microgram
µg/L	microgram per liter
mL	milliliter
ppm	parts per million
ppb	parts per billion
M	molar solution
N	normal solution

16.2 Definitions

- 16.2.1 Initial Precision and Recovery (IPR)—Four aliquots of the LRB spiked with the analytes of interest and used to establish the ability to generate acceptable precision and accuracy. An IPR is performed the first time this method is used and any time the method or instrumentation is modified.
- 16.2.2 Laboratory Control Sample (LCS)—An aliquot of LRB to which a quantity of the analyte of interest is added in the laboratory. The LCS is analyzed like a sample. Its purpose is to determine whether the methodology is in control and whether the laboratory is capable of making accurate and precise measurements.
- 16.2.3 Laboratory Reagent Blank (LRB)—An aliquot of reagent water and other blank matrix that is treated like a sample, including exposure to all glassware, equipment, and reagents that are used with other samples. The LRB is used to determine if the method analyte or other interferences are present in the laboratory environment, reagents, or apparatus.

- 16.2.4 **Matrix Spike/Matrix Spike Duplicate (MS/MSD)**—An aliquot of an environmental sample to which a quantity of the method analyte is added in the laboratory. The MS/MSD is analyzed like a sample. Its purpose is to determine whether the sample matrix contributes bias to the analytical results. The background concentration of the analyte in the sample matrix must be determined in a separate aliquot, and the measured values in the MS/MSD must be corrected for the background concentration.
- 16.2.5 **Method Detection Limit (MDL)**—The minimum concentration of a substance that can be measured and reported with 99% confidence that the analyte concentration is greater than zero.
- 16.2.6 **Minimum Level (ML)**—The level at which the entire analytical system will give a recognizable signal and acceptable calibration point, taking into account method-specific sample and injection volumes.
- 16.2.7 **Ongoing Precision and Recovery (OPR)**—See Section 16.2.2, “Laboratory Control Sample.”

17.0 Figures




Figure 1. Detailed Flow Diagram for Orthophosphate by SFA on a Flow Solution IV, Cartridge Part #A002596




Figure 2. Detailed Flow Diagram for Orthophosphate by FIA on a Flow Solution IV, Cartridge Part #A002596

Results were obtained under optimal operating conditions. Actual results may vary depending on sample introduction, cleanliness of sample containers, reagent purity, operator skill, and maintenance of instruments.

DOWFAX is a registered trademark of Dow Chemical.
Teflon is a registered trademark of E.I. du Pont de Nemours.
Flow Solution is a registered trademark of OI Analytical
WinFLOW is a trademark of OI Analytical.

Copyright 2001, OI Analytical, College Station, TX 77842.



P.O. Box 9010
College Station, Texas 77842-9010
Tel: (979) 690-1711 · FAX: (979) 690-0440