

1  
2  
3  
4  
5  
6  
7  
8  
9



# Rhithron Associates, Inc.

***Identification and enumeration of benthic diatoms using a voucher  
flora approach adapted for production taxonomy:***

***National Ecological Observatory Network (NEON)***

***DRAFT***

***Revision 3.0 2024***

***Corporate Approval***

\_\_\_\_\_  
*Wease Bollman  
President/Chief Biologist  
Rhithron Associates, Inc.  
33 Fort Missoula Road  
Missoula, MT 59804*

\_\_\_\_\_  
*Date*

\_\_\_\_\_  
*Sean Sullivan  
Operations Officer  
Rhithron Associates, Inc.  
33 Fort Missoula Road  
Missoula, MT 59804*

\_\_\_\_\_  
*Date*

Document name: Identification and enumeration of benthic diatoms using a voucher flora approach adapted for production taxonomy: National Ecological Observation Network (NEON)	Document version date: March,2024 Effective date: Date of last signature
Document number: RAI-SOP-D24v2.0	Issued by: Rhithron Associates, Inc. Chief Biologist

10

11	<b>Table of Contents</b>	
12	List of Figures.....	3
13	List of Tables.....	3
14	Executive Summary.....	4
15	Introduction .....	5
16	Diatom Project Workflow .....	5
17	Laboratory Organization .....	6
18	Health and Safety Warnings .....	7
19	Project Set-up .....	8
20	Sample intake, inventory, and chain of custody.....	8
21	Procedures .....	8
22	Sample storage and transfer to Phycology Department .....	9
23	Sample Processing and Diatom Slide Preparation.....	10
24	Materials and equipment.....	10
25	Subsampling, digestion, and slide prep .....	11
26	Sample processing QA/QC plan .....	17
27	Identification and Enumeration.....	17
28	Materials and equipment.....	18
29	Slide analysis .....	18
30	Digital images and voucher floras.....	23
31	Updates to client voucher flora .....	25
32	Metadata for image reference collections.....	26
33	Quality Assurance and Quality Control.....	27
34	QA/QC summary .....	27
35	Internal re-identification of identified transect(s) .....	28
36	External re-identification of identified transect(s) by external phycologist.....	30
37	Data Review .....	31
38	Instrument and equipment inspection, maintenance and calibration .....	32
39	Data Delivery and Sample Material Shipment.....	32

40 References ..... 34

41

42

43

44 **List of Figures**

45 **Figure 1.** Project workflow for diatom sample analysis. .... 6

46 **Figure 2:** Rhithron personnel flow chart for processing diatom projects. .... 7

47 **Figure 3 .**Project Inventory Report benchsheet example..... 12

48 **Figure 4.** Digester components. .... 13

49 **Figure 5.** Set-up for drying a low dilution and a high dilution cover slip for each sample..... 15

50 **Figure 6.** Arrangement of slides on hot plate in preparation for Naphrax addition and cover slip  
51 gluing..... 16

52 **Figure 7.** Example of sample processing metadata, entered into the Processing Data tab in  
53 Rhithron diatom data entry interface EPIC..... 17

54 **Figure 8.** Example of scribed lines on cover slip..... 19

55 **Figure 9.** Example of transect metadata entered into Sample Comments tab in Rhithron diatom  
56 data entry interface EPIC. .... 19

57 **Figure 10.** Example of a field of view with arrows marking diatoms partially in view 20

58 **Figure 11.** Labeling and notes template for a taxon in a voucher flora plate. .... 24

59 **Figure 12.** Example of an on-going client voucher template.. .... 26

60 **Figure 13.** *Example of metadata entered for a new taxon.*..... 27

61 **Figure 14.** Internal diatom QC analysis example for re-identification of same transect(s). .... 29

62 **Figure 15.** Data review flow chart. .... 32

63

64 **List of Tables**

65 **Table 1:** No aggregation of non-target taxa ..... 21

66 **Table 2:** Aggregation of non-target taxa ..... 22

67 **Table 3:** Quality control parameters and associated measurement quality objectives. .... 30

68

## 69 **Executive Summary**

70 Rhithron began adapting a voucher flora approach for identification and enumeration of  
71 diatoms starting in 2020. Motivation to regularly incorporate a voucher flora into Rhithron  
72 projects for clients was inspired by a series of publications and workshops (Bishop *et al.* 2017;  
73 Lee *et al.* 2020; Potapova *et al.* 2022; Spaulding *et al.* 2021; Tyree *et al.* 2018, 2020; USEPA  
74 2017). Detailed methods outlined how to use open-source software (GNU Image Manipulation  
75 Program, GIMP) ([www.gimp.org](http://www.gimp.org)) to organize digital images of diatom taxa into a voucher flora.  
76 Rhithron phycologists saw the enormous reference value of these vouchers for improving  
77 accuracy and precision of diatom identifications on a day-to-day basis. Rhithron subsequently  
78 adapted and modified the GIMP-based voucher flora approach to suit the unique demands of  
79 production taxonomy. Specifically, many of Rhithron’s clients have multi-year contracts that  
80 include shipments of samples scattered throughout a contract time frame. Often clients provide  
81 samples from the same locations throughout these multi-year contracts. It was essential to  
82 establish a continuous client voucher that could be referenced and updated as samples were  
83 processed over the life of the contract. Rhithron phycologists found it effective and efficient to  
84 build a client specific voucher and update it regularly as samples were analyzed. Additionally,  
85 Rhithron phycologists found it very useful to have the taxonomic nomenclature, taxonomic  
86 concepts, references, and comments noted directly in the voucher. Finally, the voucher flora  
87 approach added an additional element of quality assurance that was reflected in results from  
88 routine internal diatom QC (e.g. re-analysis of same transects by QC phycologist). Rhithron’s  
89 approach to constructing and managing an on-going voucher flora for clients is outlined in this  
90 document for consideration by the NEON working group.

91

92

93 The suggested citation for this document is:

94 Davis, C.J.; Sullivan, S.P.; Vander Meer, D.; Malik, H.; Pedraza-Garzon, E. 2024. Identification and  
95 enumeration of benthic diatoms using a voucher flora approach adapted for production  
96 taxonomy: National Ecological Observatory Network (NEON). RAI-SOP-D24V1.0. Rhithron  
97 Associates, Inc. Missoula, Montana.

98

## 99 **Introduction**

100 The Rhithron laboratory processes and identifies diatoms and soft algae in periphyton and  
101 phytoplankton samples from clients throughout North America. Data generated from these  
102 samples must be consistently and reliably generated to support client goals, which are typically  
103 related to the assessment of water quality and habitat integrity in surface water systems. The  
104 laboratory methods and protocols applied to samples vary and depend on client-specifications  
105 and project goals. Thus, samples must be handled with the utmost attention and care, and the  
106 client-specified protocol, including the required taxonomic resolution, must be faithfully  
107 followed. **This document outlines the procedures used in the analyses of diatom samples**  
108 **collected as part of the National Ecological Observatory Network (NEON) program.** These  
109 procedures are adapted from Rhithron’s Periphyton SOP ver.21.2.b(Rhithron 2021). Rhithron’s  
110 internal data quality objectives (see QAP Ver 21.1.a) are maintained by following these  
111 procedures.

112


113 All personnel working on the NEON project samples, in any capacity, must review the version-  
114 controlled Standard Operating Procedure (SOP) in its entirety prior to working on the project.

115 Additionally, if more than two months lapse without a review of the NEON SOP then each  
116 person shall review the SOP in its entirety.

## 117 **Diatom Project Workflow**

118 Samples received at Rhithron enter an organized workflow specifically for a diatom analysis  
119 project (Figure 1). **A project is defined as a shipment of client samples received by Rhithron**  
120 **(e.g., 15-25 samples) that will proceed through the workflow.** Clients typically have several  
121 Rhithron projects completed annually to accommodate their field sampling activities. The  
122 workflow diagram in Figure 1 illustrates the timing of key activities for each client’s project.  
123 Each activity is presented in detail within this SOP.

124




125  
126 **Figure 1.** *Project workflow for diatom sample analysis.*

127

## 128 **Laboratory Organization**

129 The organizational chart in **Error! Reference source not found.** shows the Rhithron personnel  
130 responsible for the various tasks associated with the workflow and illustrates the pathways of  
131 communication that assure the quality of Rhithron’s work.

132



133

134 **Figure 2:** Rhithron personnel flow chart for processing diatom projects.

## 135 Health and Safety Warnings

136 Health and safety procedures must be followed to prevent laboratory injuries. All proper  
137 personal protection clothing and equipment (e.g., lab coat/apron, protective eyewear/goggles,  
138 nitrile/rubber gloves, fume hood) must be worn or applied. When working with potentially  
139 hazardous chemicals (e.g., glutaraldehyde) or biological agents (benthic organisms or sediment)  
140 avoid inhalation, skin contact, eye contact, or ingestion. If skin contact occurs remove clothing  
141 immediately and wash/rinse thoroughly. Wash the affected skin areas thoroughly with large  
142 amounts of soap and water. All staff performing duties associated with chemicals in this SOP  
143 shall review the Rhithron Health and Safety Manual (Rhithron 2022) and review all Safety Data  
144 Sheets for any of the chemicals they will handle. Safety Data Sheets are maintained in the

145 Rhithron laboratory Safety Plan.

146

147

## 148 **Project Set-up**

### 149 ***Sample intake, inventory, and chain of custody***

150 Rhithron is notified of sample shipment through the CLANotify system within the NEON  
151 framework, and receives the manifest and receipt forms at the time of shipping notification.  
152 Upon delivery the Intake Specialist takes custody of the delivered samples. All samples are  
153 shipped directly to the Rhithron laboratory, no off-site location is used.

154

155 The Intake Specialist assesses the condition of the samples and preservative needs, makes sure  
156 that all samples correspond with the NEON shipment manifest, and that all expected parts of  
157 the delivery have arrived safely. The Intake Specialist is also responsible for the completion of  
158 sample intake procedures, including generation of the project inventory report, chain of  
159 custody records, NEON receipt form and transfer of custody to the Phycology Department upon  
160 completion of the intake and inventory procedures.

161

### 162 ***Procedures***

163 When a project arrives, all shipping containers are opened by the Intake Specialist and the  
164 manifest is located on the Rhithron server. The manifest is referred to during all of the  
165 following steps. Each sample jar is removed, and the level and integrity of preservative is  
166 checked, and recharged with glutaraldehyde (1% final concentration volume/volume) when  
167 required (e.g. decaying matter smell detected). The label of each sample container is checked  
168 against the manifest, and each sample is marked off on the manifest as it is identified. The  
169 NEON-provided barcode label is scanned. For samples with multiple bottles, all bottles are  
170 organized together and scanned. Any leakage or damage, and any discrepancies between  
171 sample labeling and the manifest document are noted. This information is reported to the  
172 Operations Officer immediately. The Operations Officer reports this information to NEON by  
173 replying (Reply All) to the shipping email notification to immediately notify pertinent parties to  
174 the problem and also by reporting this on the sample receipt forms. Discrepancies are rectified  
175 before manifests are signed, copied, and returned to the client.

176 Once all samples have been checked against the manifest, and all discrepancies have been  
177 rectified by the client, the NEON receipt form is completed by the Intake Specialist. Upon  
178 completion, the receipt form is reviewed by the Data Specialist and uploaded to the NEON data  
179 repository.

180 The Intake Specialist creates the Rhithron inventory file (RAI Inventory) from the client  
181 manifest, which is an Excel file saved to the Rhithron server in the NEON client folder. Entries in  
182 this file include the site name, client sample identifiers, replicate numbers, sample collection  
183 dates, the number of bottles in each sample (e.g., 1 of 2, 2 of 2, etc.), and any other  
184 distinguishing data. The RAI Inventory contains a number of worksheets: 1) the “*Client*”  
185 worksheet contains client Chain of Custody (COC) information, and any other information  
186 gathered by the Intake Specialist from bottle labels, 2) the “*Pre-check-in*” worksheet contains  
187 client information that is copied into the Rhithron Associates, Inc. Laboratory Information  
188 System (RAILIS) format for subsequent upload, and the creation of sample labels and RAI  
189 numbers (internal identifiers) assigned to each sample, 3) the “*Check-in*” worksheet contains  
190 client-corrected information consistent with the bottle labels, and 4) the “*Sample Upload*”  
191 sheet is used by the Intake Specialist to upload sample information to Rhithron’s proprietary  
192 Electronic Periphyton and Invertebrate Cloud software (EPIC). The Intake Specialist files the  
193 physical copy of the client’s COC document.

194 The Intake Specialist uploads the internal inventory into RAILIS and prints the sample labels  
195 output, which is generated by RAILIS. The sample label output consists of sample identifier  
196 labels for each of the sample containers, and labels for the processed diatom slides. Each label  
197 type includes the unique internal identifiers (RAI numbers) assigned to each sample. The Intake  
198 Specialist attaches the appropriate sample identifier labels to the sample bottles during check-  
199 in.

200 The Intake Specialist takes samples to the storage site. The location of samples within a project  
201 are recorded as coordinates within the storage site (e.g., Fridge 1: Shelf A). Location  
202 coordinates and preservative type are recorded in RAILIS. RAILIS generates a Project Inventory  
203 Report benchsheet and sample labels for use during processing; these are printed when the  
204 samples are about to be processed. The Intake Specialist notifies the Phycology Department of  
205 the status of the samples and their location in the storage site.

### 206 ***Sample storage and transfer to Phycology Department***

207 Periphyton samples awaiting processing are stored under refrigeration in the processing  
208 laboratory storage site. RAILIS enables staff to track custody transfers, current sample location,  
209 and the progress of a project through laboratory analysis. Once sample intake procedures are  
210 complete, custody of the project and samples is transferred to the Phycology Department  
211 which is responsible for the security of samples stored while awaiting processing.

212  
213 Project inventory reports are delivered to a current projects folder in the Phycology  
214 Department office area. All phycologists are notified of the project sample locations by the

215 Intake Specialist via email. Stored samples and sample fractions are checked monthly by  
216 phycologists to ensure sample integrity, and preservatives are added as needed.

217

## 218 **Sample Processing and Diatom Slide Preparation**

219 Goals of diatom sample processing includes digestion, which produces a “cleaned”  
220 subsample portion, free of organic matter, and containing a diatom assemblage representative  
221 of the sample. The subsample is then used to prepare a slide mount with an even distribution  
222 of this assemblage. Samples are processed in a standardized, repeatable manner. In general,  
223 the digestion of diatom sample material follows standard methods outlined in the Academy of  
224 Natural Sciences SOP for the National Assessment of Water Quality (NAWQA) program (Charles  
225 *et al.* 2002). Diatom processing procedures are applied to achieve the following outcomes:

- 226 • The client-specified protocol is followed precisely. Client specifications are reviewed by the  
227 Operations Officer, phycologists, and Intake Specialist before the processing of samples. If  
228 there is any question or uncertainty about any procedure or protocol detail, the Operations  
229 Officer is consulted for clarification before proceeding.
- 230 • The appropriate paperwork (i.e., sample labels, slide labels, Project Inventory Report benchsheet)  
231 is associated with the correct project samples.
- 232 • All data pertinent to the processing procedure, including initial sample volume,  
233 processed volume, and final digested volume are recorded and entered into the database  
234 appropriately.
- 235 • Cross-contamination among samples does not occur.

### 236 **Materials and equipment**

237 50 mL centrifuge tubes

238 RAILIS-generated sample labels

239 Liquinox® soap

240 Microwave digestion system (Milestone ETHOS EZ), 10 rotor digestion vials

241 70% nitric acid

242 Bottle-top dispenser (BRAND® Dispensette® S, Organic, Analog-adjustable, 2.5-25 mL volume  
243 range)

244 Distilled water (using Waterwise 9000 Countertop Distiller), in wash bottle

245 Graduated cylinder

246 Ultrasonic bath (Garosa 1.3L Model 113HTD)

247 Project Inventory Report benchsheet (pre-printed and specific to project)

248 Hot plate, temperature controlled

249 Glass microscope slides (2.5x7.5cm)

250 Glass cover slips (round 18mm diameter), No.1 thickness

251 Transfer pipettes, disposable (7mL capacity, 0.5mL graduation)

252 Razor blades, single edged

253 Naphrax™ mounting medium

254 Forceps  
255 Drying pan  
256 20 mL glass scintillation vials  
257 90% acetone (PHARMCO, general use HPLC-UV grade), in 50 mL glass vial  
258 Compound microscope (Leica DM1000), equipped with 10X and 63X objectives  
259 PPE (gloves, goggles, lab coat/apron)  
260 Fume hood  
261 pH indicator strips

262

### 263 ***Subsampling, digestion, and slide prep***

264

265 **Step 1:** Label 50 mL centrifuge tubes with the preprinted RAILIS-generated label for each  
266 diatom sample and place them in a centrifuge tube rack. Organize project samples into batches  
267 of 10, which is the capacity of the digester. For small projects, create batches that divide the  
268 samples as close to evenly as possible. For example, a project of 13 total samples will be a batch  
269 of 6 and a batch of 7.

270 **Step 2:** Homogenize periphyton sample in its container(s) by thoroughly agitating by hand for  
271 one minute to dislodge diatom epiphytes from filamentous algae and to mix sample. Pour the  
272 entire periphyton sample into an appropriately sized graduated cylinder to obtain the initial total  
273 volume. Record the initial sample volume on the Project Inventory Report benchsheet in “Initial  
274 vol (mL)” column (**Error! Reference source not found.**).

275 **Step 3:** Pour a subsample of 25 mL (minimum 5 mL) into the corresponding labeled centrifuge  
276 tube. Note and enter the volume on the Project Inventory Report benchsheet in “Processed vol  
277 (mL)” column (**Error! Reference source not found.**).

278 **Step 4:** Between samples, ensure that there is no cross-contamination by thoroughly rinsing  
279 the graduated cylinder. Completely fill and rinse the graduated cylinder at least 3 times with  
280 distilled water.

281 **Step 5:** Record processing date (month/day/year) on the Project Inventory Sheet in “Processed  
282 Date” column (**Error! Reference source not found.**). Record any notable features of the  
283 unprocessed sample material (e.g. high detritus, sediments, green algae clump) in the  
284 “Processing Notes” column (**Error! Reference source not found.**).

285 **Step 6:** Repeat steps 2 through Step 5 until all samples have been subsampled for diatom  
286 processing.

287

Project Inventory Report				Project ID:	WADOE23BIP			
				Project Name:	Biological Integrity Mon Peri 2023			
<i>Customer Information</i>		<i>Project Management Information</i>		<i>Project Location Information</i>				
WADOE		Date Arrived: 10/23/2023		Country: USA				
Chad Larson		Date Due: 6/1/2024		State: Washington				
(360) 407-7456		No. Samples: 14		County:				
clar461@ECY.WA.GOV		No. Jars: 14		City:				
Sample Information								
RAI Number	Client ID	Date Coll.	No. Jars	Processed Date	Initial vol (mL)	Processed vol (mL)	Rotor vial #	Processing Notes
WADOE23BIP001	BIO06600-BURP15-DCE-2023-0718-10:30	7/18/2023	1	01/01/2024	980	20	1	High organics

288

289 **Figure 3** .Project Inventory Report benchsheet example, illustrating key sample information (  
 290 i.e., initial volume, processed volume, processing notes) to be recorded at the bench during the  
 291 subsampling procedure.

292 **Step 7:** Top off the centrifuge tubes to 50mL with distilled water. Spin the centrifuge tubes in the  
 293 centrifuge for 15 minutes at 2500 RPM. After centrifuging, carefully remove the supernatant from  
 294 each tube, using the vacuum spigot, and being sure to **stop vacuuming at 7.5mL in the tube to**  
 295 **avoid removing the settled material.**

296 **Step 8:** Transfer the remaining material, including the settled material, from the centrifuge tube  
 297 to a rotor digestion vial from the microwave digester. The rotor digestion vials are already labeled  
 298 from 1 to 10 (Figure 4). Using distilled water, rinse the centrifuge tube with 7-10mL, adding this  
 299 rinsate to the digestion vial. Try to obtain a total volume of approximately 15 mL of sample  
 300 material: this will affect the final acid-to-sample-volume ratio: see Step 9. Be careful not to  
 301 overdo the rinsing. Insert each rotor digestion vial into the rotor and record the position of each  
 302 sample on the Project Inventory Report benchsheet in the “Rotor vial #” column (**Error!**  
 303 **Reference source not found.**).

304 **Step 9:** For this step, wear protective equipment, including acid-resistant gloves, a lab jacket or  
 305 apron, and eye protection. This step must be performed inside the fume hood, with the sash  
 306 down. Confirm the bottle-top dispenser is set to 15mL and carefully add 15mL of 70% nitric  
 307 acid to each of the rotor digestion vials using the bottle top dispenser. The final acid-to-sample-  
 308 volume ratio should be as close to 1 (equal parts acid and sample) as possible.

309 **Step 10:** Secure the rotor digester vials into the digester blocks using a torque wrench set to 22.5  
 310 Nm.

311 **Step 11:** Place blocks in microwave digester verifying rotor vial, block, and turret position (Figure  
 312 4). Carefully secure temperature probe in block #1 (refer to Milestone ETHOS EZ manual for  
 313 proper probe installation and care). Temperature probe should fit in smoothly (i.e. does not  
 314 require force), if not, inspect probe for damage and note alignment of digestion components.



315

316 **Figure 4.** *Digester components labeled 1-10 to track sample locations during digestion process.*  
317 *Red circles indicate example of rotor digester vial #7 location in digester block #7 and placed in*  
318 *position #7 within the digester turret.*

319 **Step 12:** Run the program for the digester. The temperature will increase from room  
320 temperature to 140°C for the first 20 minutes. Over the next 40 minutes the temperature will  
321 gradually increase to 170°C. The final phase will be a 10-minute cool down period. Total run  
322 time will take 70 minutes.

323 **Step 13:** Remove the blocks individually from the digester, place them in the fume hood, and  
324 allow them to cool for 45 to 60 minutes.

325 **Step 14:** Remove the rotor vials from the blocks and transfer the digested material to the labeled  
326 centrifuge tubes using distilled water to rinse the rotor vials.

327 **Step 15:** Carefully dilute the nitric acid by topping off the centrifuge tubes with distilled water  
328 using a wash bottle. Avoid splashing and overspray. Centrifuge the samples for 15 minutes at  
329 2500 RPM. Carefully siphon the liquid down to 7.5 mL using the vacuum spigot. Repeat dilution-  
330 centrifuge-siphon steps at least 6 times, then check with pH indicator strips to ensure the acid is  
331 sufficiently neutralized to greater than 5 pH.

332 **Step 16:** Record the final volume after digestion on the Project Inventory Report benchsheet in  
333 the “Processing Notes” column (**Error! Reference source not found.**). Typically, digested diatom  
334 samples should be 7.5 mL at this step.

335 **Step 17:** Steps 17 through 27 will result in the assembly of 2 diatom slide mounts for each  
336 sample, a lower-density slide mount and a higher-density slide mount. This helps to ensure that  
337 the phycologist has a slide that can be easily analyzed. Cover slips are cleaned in 90% acetone  
338 and wiped dry with laboratory wipes prior to use. Discard any cover slips that have debris that  
339 cannot be removed. Level the hot plate. **Label a slide with a RAILIS-generated sample label AND**  
340 **write the rotor digester vial number above the sample RAI number.** Place two (2) cleaned  
341 round 18mm cover slips for each sample on the labeled slide (Figure 5). Record the rotor


342 digester vial number on the Project Inventory Report benchsheet and ***ensure that the rotor***  
343 ***digester vial number matches the RAI sample number.***

344 **Step 18:** Using a clean transfer pipette for each sample, mix the solid and liquid portions  
345 remaining in the centrifuge tubes and transfer this material to a 20mL vial labeled with the  
346 RALLIS-generated sample label. Gently sonicate clean diatom material to break up clumps by  
347 smoothly dipping vial in/out of sonication bath 10 times. Sonicator water temperature should  
348 be room temperature so temperature control should be off.

349 **Step 19:** In steps 19 and 20, the low dilution (higher-density) and high dilution (lower-density)  
350 subsamples are prepared from the material in the 20 mL vial. Visually evaluate the biomass of  
351 diatom material in the vial by holding it up to a light source, and swirling the vial to suspend the  
352 solid materials. Diatom frustules will appear as “white sparkles” while presence of sand  
353 particles will add a brown/darkening of the sample. Based on observations and experience,  
354 determine whether the diatom biomass is relatively high: water may be added to high biomass  
355 samples. If water is added, the total sample volume should not be allowed to exceed 20 mL.  
356 Update the final digested volume recorded in the Project Inventory Report benchsheet in the  
357 “Processing Notes” column (Figure 3) to reflect the addition of water. Needing to concentrate  
358 samples in cases of very low biomass is extremely rare, however centrifuging-siphoning to a 1  
359 mL can be achieved.

360 **Step 20:** Choose a low (e.g., 50) and a high (e.g., 100) dilution factor, based on observation and  
361 experience. Prepare the low dilution cover slip first, followed by the high dilution cover slip.  
362 Agitate vial thoroughly before subsampling. Mix dilutions in a clean graduated glass beaker by  
363 using a clean transfer pipette to add 1 mL of cleaned material and bring up to the final volume  
364 with distilled water. For example, 1 mL sample + 49 mL water will result in a dilution factor of  
365 50. **Confirm the location of the labeled slide that matches the diluted sample** and transfer  
366 approximately 1 mL of a diluted mixture onto a clean cover slip. **Low dilution goes on the left-**  
367 **side cover clip and high dilutions go on the right-side cover slip** (Figure 5). Rinse the graduated  
368 beaker thoroughly with distilled water between samples. Repeat dilution steps for high dilution  
369 cover slip, transferring the second diluted mixture to the right-side cover slip. In the Processing  
370 Data tab in EPIC (Figure 7), record the dilution/concentration factor (DCF) for the low dilution  
371 cover slip in the Slide Dilution Concentration Factor field, and record the DCF for the high  
372 dilution cover slip in the Notes field.

373 **Step 21:** Left-over cleaned material in the 20 mL vial is stored at room temperature in a project-  
374 labeled box. Lids are left loose to allow water to evaporate leaving dried diatom frustules.  
375 Cleaned material in barcoded vials are checked after 10 days to determine if additional time is  
376 needed for desiccation. The lids of dried materials are secured at this time.



377

378 **Figure 5.** Set-up for drying a low dilution and a high dilution cover slip for each sample. In this  
379 example, 2 samples are shown (058, 059) for this client's project (ECO A22NL05).

380 **Step 22:** Allow the cover slips to dry overnight at room temperature (~70°F). Examine the cover  
381 slips the following morning to see if all liquid has apparently evaporated. Temperature and  
382 humidity will vary by season/day, and this can impact drying times.

383 **Step 23:** Scan dried cover slips under the microscope to evaluate valve density and dispersion,  
384 and the quality of the dried subsample. At 10X magnification, assess for the presence of  
385 sediment, clumping of frustules, residue, or other problems. Using the 63X objective, evaluate  
386 valve density. Ideally, each field of view (FOV) at 63X should have approximately 20 diatom  
387 valves visible. Additional processing steps, including additional sample rinses to remove  
388 residues or dilutions for high valve density or excessive sediment, may be needed.

389 **Step 24:** With a hot plate set to the lowest setting (i.e. warm to the touch), gently heat the  
390 dried cover slips for 1 hour to remove any remaining moisture, which may not be visible to the  
391 naked eye. Using the Project Inventory Report benchsheet, label slides on the left side with the  
392 corresponding RAILIS-generated sample label and **write the dilution factor above RAI numbers.**

393 **Carefully refer to the benchsheet and to the previously-labeled slide to ensure that the RAI**  
394 **numbers on all of the slides corresponding to the sample are in agreement.** Place labeled  
395 slides on a hot plate set to 250°C inside the fume hood, arranging along the edge in the  
396 following order: project sample 001 low dilution, project sample 001 high dilution, project  
397 sample 002 low dilution, project sample 002 high dilution, etc. (Figure 6). Using a transfer  
398 pipette, place a small pea-size drop of Naphrax on each heated slide. Naphrax should be easily  
399 pipetted, having a viscosity similar to glycerin. If the Naphrax is too viscous, it is likely that the  
400 toluene has evaporated, and new Naphrax is needed.

401 **Step 25:** Carefully invert the cover slip using forceps so that the dried diatom material is placed  
402 on the heated slide, on top of the drop of Naphrax.



403

404 **Figure 6.** Arrangement of slides on hot plate in preparation for Naphrax addition and cover slip  
405 gluing. In this example, 2 project samples are shown (058, 059) for this client’s project  
406 (ECO A22NL05), each having slides with a low dilution (50 DCF) and high dilution (100 DCF).

407 **Step 26:** Allow the slides and cover slips to sit on the hot plate briefly (typically 3-5 minutes) to  
408 allow the toluene to boil out of the Naphrax. Bubble formation should cease when all toluene  
409 has been removed. Remove the slides from the heat and gently press the cover slips with forceps  
410 to remove any bubbles in the Naphrax.

411 **Step 27:** Allow the slides to cool for a few minutes and carefully scrape excess Naphrax from  
412 around the cover slips using a razor. The Naphrax should be brittle and easily removed. If  
413 otherwise (e.g., sticky, soft), return the slide to the hot plate for additional heating to remove  
414 toluene.

415 **Step 28:** Scrub all labware using Liquinox © solution and rinse with tap water followed by triple  
416 rinse with distilled water.

417 **Step 29:** Enter metadata from sample processing into the Processing Data tab (Figure 7) in EPIC,  
418 the diatom data entry interface. Enter the initials of the processor, the processing date, initial  
419 volume, processed volume, slide condition, DCF value of “best” slide and general processing  
420 notes (e.g., additional slide made 50 DCF).

421

Diatom Taxonomist	
Processing Tech:	CD
Processing Date:	02/14/2024
Initial Volume:	50
Processed Volume:	25
Slide Condition:	High sediment
Slide Dilution Concentration Factor:	25

**Processing Notes**

Additional slide 50 DCF

422  
423 **Figure 7.** Example of sample processing metadata, entered into the Processing Data tab in  
424 Rhithron diatom data entry interface EPIC.

425

### 426 **Sample processing QA/QC plan**

427

428 Sample volumes and acid volumes are carefully monitored by the processing staff in order  
429 to maintain equivalent processing procedures for all samples. Great care is taken to ensure  
430 that contamination does not occur during processing by thoroughly cleaning all  
431 equipment, and using the microwave digester, which isolates each sample during the  
432 digestion procedure.

433

### 434 **Identification and Enumeration**

435 All analysis of diatom samples is performed by Rhithron staff phycologists. Phycologists have 5  
436 to 20 years of experience in the species identification of diatoms, and each holds Society for  
437 Freshwater Science (SFS) Level 1 diatom genus or level 2 diatom species taxonomic  
438 certifications. The following methods are used to ensure reproducible diatom identification  
439 and enumeration data from processed sample material. Phycologists analyzing the diatom  
440 slides should note this very important project data element:

- 441 • New clients having their first project voucher flora completed will have that project voucher  
442 saved as the on-going client voucher that will be referenced, and added to, in subsequent  
443 projects. The first project voucher is a time-consuming process as all taxa encountered are

444 effectively “new” and documented following the standardized approaches outlined in this  
445 Identification and Enumeration section.

446 • **Identified taxa are assigned species names at the bench by referencing images in a client**  
447 **voucher flora. If a specimen does not match any taxon in the client voucher, the phycologist**  
448 **captures an image and assigns a name, as described in Step 5 below.** Once the current project  
449 is completed, all new images and taxa are compiled into the project voucher flora and reviewed  
450 by at least two Rhithron phycologists. Phycologists discuss disagreements and review literature  
451 as a group. A consensus is generally reached, with an occasional “majority rule” approach, and  
452 on rare instances a “tied” result requiring consultation with an external phycologist for input.  
453 At this point, initial bench names are subject to change. **New taxa in the project voucher flora**  
454 **are added to the on-going client voucher for reference when identifying diatoms in future**  
455 **client projects.**

456

### 457 ***Materials and equipment***

458 Compound Microscope (Leica DM2500), equipped with differential interference contrast (DIC)

459 Diamond scribe, mounted on objective turret

460 100X oil objective (1.4NA)

461 10X, 20X,40X objectives

462 Oil condenser

463 Digital camera (Leica DMC2900)

464 Networked computer located at the microscope for access to EPIC data entry software

465 GNU Image Manipulation Program (GIMP)

466 Measurement calibration slide (1mm ruler, 0.01mm resolution)

467

### 468 ***Slide analysis***

469 **Step 1:** For a sample, select one of the prepared slides (high density or low density) for analysis.

470 Use a Leica DM 2500 compound microscope equipped with DIC, an oil condenser, and 100X

471 objective (1.40NA, 1000x magnification) for identification and enumeration.

472 **Step 2:** Scan the slide using coarse (100X) magnification to confirm adequate density and valve

473 dispersion. Also note areas of optical distortions that should be avoided, such as areas near the

474 cover slip edges. Additionally, the second slide can also be scanned as an alternative with a

475 higher/lower DCF value. Confirm that the correct DCF value is entered in EPIC for the chosen

476 slide (Figure 7). Sometimes neither slide is deemed adequate, and the phycologist can choose

477 to re-process the sample, obtaining a different DCF, in order to obtain a more readable diatom

478 slide.

479 **Step 3:** Using the diamond scribe, scribe a vertical transect line in an area free of optical

480 distortions and avoids the cover slip edges. A horizontal “start” line is scribed to mark the

481 beginning of the vertical transect and the count proceeds along the right side of the transect

482 line, as viewed through the microscope (Figure 8). Another horizontal line is scribed at the “end”

483 of the transect when counting along the transect ceases. Record transect location, start, and end

484 locations using the X/Y coordinates on the stage Vernier scale (e.g., transect 1 start 38/109,  
485 transect 1 end 38/93). Several transects may be required to attain the target count of 600  
486 valves; each transect should be scribed as described above. Record these data in the Sample  
487 Comments tab, in the Comments field within EPIC (Figure 9). Record the total length and  
488 number of transects when the sample identification and enumeration are completed (Figure 9).  
489  
490



491  
492 **Figure 8.** Example of scribed lines on cover slip. Vertical scribed lines are viewed top-to-bottom  
493 on the right-hand side of the transect line. Starting and ending points are indicated by  
494 horizontal scribed lines. RAI slide label on the left side of coverslip and NEON unique barcode to  
495 the right of the coverslip.


496

497  
498 **Figure 9.** Example of transect metadata entered into Sample Comments tab in Rhithron diatom  
499 data entry interface EPIC. Fields of view are not recorded using the transect method.  
500 Dilution/Concentration Factor is not recorded here.

501 **Step 4:** A diatom valve or frustule (2 valves) is counted if greater than half of the specimen is  
502 within the FOV (Figure 10); this includes specimens lying at the edge of the FOV as well as  
503 broken valves. A fragment is given a count of 1 if over half of the valve is visible within the FOV,  
504 including the central area. For raphid species, both central raphe endings should be visible.

505 **Step 5: Identify diatoms to the lowest possible taxonomic level, generally to species, by**

506 **referencing the client voucher flora. As identifications proceed, enter taxon names, counts,**  
507 **and comments into the Taxon, Count, and Comments fields, respectively, in the Taxonomic**  
508 **Data Entry interface, which is part of the EPIC software. The Taxonomic Data Entry interface**  
509 **is error-protected for taxonomic nomenclature. For specific EPIC data entry procedures, see**  
510 **the EPIC Instruction Manual for Rhithron Phycologists**(Rhithron 2014). Taxa that are not in the  
511 on-going client voucher flora are documented with images (see recommendations in section  
512 **Digital images and voucher floras**), identified following standard taxonomic references, and  
513 assigned current nomenclature (e.g., Spaulding et al 2021, Cantonati et al. 2017). Rhithron’s  
514 taxonomic resource library includes an extensive collection of books and documents in hard  
515 copy/electronic form on the Rhithron network server. Taxa that are not found in this resource  
516 library are identified to genus and assigned provisional names (e.g., *Navicula* sp1) in the  
517 Comments field for the taxon entry. Occasionally, genus-level identification is uncertain (e.g.,  
518 when the specimen is obscured); in these cases the specimen is designated as “undetermined”,  
519 and names such as “Undetermined Centric” and “Undetermined Pennate” are used. These  
520 names should be used sparingly, and, ideally, only provisionally. Phycologists should use a  
521 combination of sample images, additional scanning of slides, and experience with taxa and  
522 existing vouchers to help achieve target-level identifications for these troublesome specimens.  
523  
524



525  
526 **Figure 10.** *Example of a field of view with arrows marking diatoms partially in view.*  
527 *Specimens marked with a solid black arrow would be included in the count, whereas those*  
528 *with a white arrow would be excluded.*

529  
530 **NEON WORKING GROUP NEEDS TO MAKE A DECISION ON HOW TO HANDLE NON-TARGET**  
531 **TAXA ENUMERATION ( we suggest either going with step 6a OR step 6b)**  
532 **Step 6a: (NO AGGREGATION)** Diatom valves that are oriented in girdle (side) view, are damaged,  
533 or are obscured because of tilting or debris may not be able to be identified to target-level.  
534 These must be documented consistently during the analysis. Phycologists can make a temporary

535 genus-level entry (e.g., *Nitzschia*) with a descriptive note (e.g., girdle, 20x3um cf *N. fonticola*) in  
 536 the Comments field (See Table 1) in EPIC. Ideally, the identity of such specimens can be resolved  
 537 as more specimens are confidently identified during the analysis. Situations do arise where  
 538 assigning these girdle/obscured/damaged valves to a species cannot be made with confidence. A  
 539 common example would be *Achnantheidium* girdle views where multiple species from the genus  
 540 can have very similar girdles. For such a case, enter “*Achnantheidium*” and make a note “girdle” in  
 541 the Comments field (Table 1). All girdle views that look similar will be counted in this entry  
 542 throughout the entire analysis of the slide (e.g., 300 girdle valves counted). *Rhoicosphenia*,  
 543 *Gomphonema*, and *Eunotia* girdle views, *Planothidium* and *Psammothidium* rapheless valves,  
 544 and *Cocconeis* raphe valves are all possible situations that may require this approach. **Decisions**  
 545 **on how to analyze non-target taxon entries (e.g., aggregations) will be left to the end user of**  
 546 **the dataset.**

547  
 548 **Table 1:** *No aggregation of non-target taxa (i.e., girdle views not added proportionally to species*  
 549 *entries)*

Taxon	Count	Comments
<i>Achnantheidium minutissimum</i>	90	
<i>Achnantheidium rivulare</i>	10	
<i>Achnantheidium</i>	300	Girdle

550  
 551 **Step 6b: (AGGREGATION)** Diatom valves that are oriented in girdle (side) view, are damaged, or  
 552 are obscured because of tilting or debris may not be able to be identified to target-level. These  
 553 must be documented consistently during the analysis. Phycologists can make a temporary  
 554 genus-level entry (e.g., *Nitzschia*) with a descriptive note (e.g., girdle, 20x3um cf *N. fonticola*) in  
 555 the Comments field for the taxon in EPIC. Ideally, the identity can be resolved as more species  
 556 are confidently identified during the analysis. Situations do arise where assigning these  
 557 girdle/obscured/damaged valves to a species cannot be made with confidence. A common  
 558 example would be *Achnantheidium* girdle views where multiple species from the genus can have  
 559 very similar girdles. For such a case, enter a temporary entry of “*Achnantheidium*” and make a  
 560 note “girdle” in the Comments field for that taxon. All girdle views that look similar will be  
 561 counted in this entry throughout the entire analysis of the slide (e.g., 300 counted). In this  
 562 example, two *Achnantheidium* species are identified in the sample (e.g., *A. minutissimum* and  
 563 *A. rivulare*): the *Achnantheidium* in girdle view could be either of these species. Upon completion  
 564 of sample identification to the specified target count, genus-level girdle view counts are assigned  
 565 to the 2 species identified, in proportion to the relative abundance for each confirmed species.  
 566 For example, if *A. minutissimum* and *A. rivulare* had counts of 90 and 10 respectively, and if there  
 567 were 50 genus-level counts of *Achnantheidium* “girdle”, 45 of these 50 counts are assigned to *A.*  
 568 *minutissimum*, and the remaining 5 are assigned to *A. rivulare*. *A. minutissimum* now has a

569 count of 135, and *A. rivulare* has a count of 15, and the genus-level *Achnantheidium* “girdle” entry  
570 is deleted (Table 2). *Rhoicosphenia*, *Gomphonema* and *Eunotia* girdle views, *Planothidium* and  
571 *Psammothidium* rapheless valves, and *Cocconeis* raphe valves are all possible situations that may  
572 require this outlined approach. High numbers of girdle/obscured counts that are processed this  
573 way should be noted in the taxon entry Comment field (e.g., added 45 *Achnantheidium* girdles,  
574 Table 2).

575  
576

577 **Table 2:** *Aggregation of non-target taxa (i.e., girdle views added proportionally to species*  
578 *entries)*

Taxon	Count	Comments
<i>Achnantheidium minutissimum</i>	135	Added 45 <i>Achnantheidium</i> girdles
<i>Achnantheidium rivulare</i>	45	Added 5 <i>Achnantheidium</i> girdles

579

580 **Step 7:** Count each specimen along the scribed transect(s) until the target count is achieved for  
581 each sample. Enter taxon names, counts, and other information as you identify and count  
582 specimens. Continue identifications and counts until any one of the following criteria is met:  
583 1) the target count is reached OR 2) a maximum of **5 transects is analyzed, OR 3) 4 hours,**  
584 **excluding time spent imaging and/or learning new taxa, is spent in the analysis.**

585 **Step 8:** When counts and identifications are complete, attach the barcode sticker to the  
586 analyzed slide, on the right side of the coverslip (opposite Rhithron sample label) (Figure 8).  
587 Barcode numbers for each sample are recorded in the shipping manifest by the Intake Specialist  
588 when slides are sent to the Biorepository at Arizona State University. Additional slides that  
589 were not analyzed will be archived at the NEON biorepository at ASU.

590  
591  
592  
593  
594  
595  
596  
597  
598  
599  
600  
601  
602  
603  
604  
605

606 ***Digital images and voucher floras***

607 Images should be taken during analysis to document new taxa and improve identification  
608 accuracy and precision. In general, the following recommendations should be considered:

- 609
- 610 • Images are needed for **all new taxa that are not documented in the client diatom voucher**. If  
611 a diatom voucher flora does not yet exist for a client, the first project voucher created is  
612 copied and saved as the client voucher that will be referenced and added to for future  
613 projects. A subfolder for the new client voucher flora should be created in the appropriate  
614 client folder on the server. Use a descriptive name for the subfolder, such as NEON23\_Diatom  
615 Voucher Flora.
  - 616 • Record the cover slip XY coordinates for each new taxon in the Reference Collection tab (see  
617 Figure 10). For more details, see the section: Metadata for image reference collections,  
618 below.
  - 619 • In addition to new taxa, images should be taken for 1) difficult taxonomic groups (e.g.,  
620 *Fragilaria*, *Gomphonema*, *Nitzschia*, *Sellaphora*), 2) any taxa that require further scrutiny and  
621 “keying out” in references, or 3) specimens that closely fit existing species concepts in the  
622 client voucher but have attributes that “expand” the concept.
  - 623 • Images should be of sufficient quality to illustrate the important morphological characteristics  
624 that define the species. Ask yourself, “Could I really ID this specimen to species with this  
625 image?”. If in doubt, share the image with another phycologist to get some feedback. Bad  
626 images are not worth reviewing or adding to a voucher.
  - 627 • Capture images using the grey scale setting in the camera software.
  - 628 • Always include a scale bar within the image. Place the scale bar in the lower right corner of  
629 the image. The scale bar should be 10 microns long and vertically oriented.
  - 630 • Name the image file in the following format Taxon\_Sample ID\_Magnification\_YearMonthDay.  
631 For example: *Gomphonema parvulum*\_NEON24001\_1000x\_20240721.
  - 632 • Export images into the project subfolder (e.g., NEON24\_Raw Diatom Images) on the Rhithron  
633 server.

634


635 **Step 1:** Review raw diatom images for the project in the project subfolder on the Rhithron  
636 server. Decide which images/taxa should be processed and included in the project voucher  
637 flora. Priority is given to taxa that are 1) not currently documented in the on-going client  
638 voucher flora or 2) abundant taxa in extremely difficult taxonomic groups. Inclusion in the  
639 project voucher flora is not limited to these taxa; phycologists can include as many taxa as they  
640 like in the project voucher. Remember to occasionally check the taxonomic concepts in the on-  
641 going client voucher flora to make sure the species concept has not drifted during current  
642 analysis.

643 **Step 2:** Process selected images and arrange them into the project voucher flora using the  
644 GIMP software. Details for how to process images and create a voucher are outlined in  
645 materials provided by DONA Diatom Web Academy workshop, May 2020  
646 (<https://diatoms.org/news/online-voucher-flora-workshop>)(Lee *et al.* 2020). Arrange taxa  
647 images by genus or morphology (depending on number of images) into separate flora plates  
648 (e.g., NEON24\_Achnanthidium\_plate1.psd, NEON24\_Surirelloids\_plate1.psd ). Differences in  
649 camera age and quality on Rhithron’s scopes requires that GIMP-processed images (Scaled  
650 Image in X and Y resolution in GIMP) should have a standardized resolution of 72 dpi.

651 **Step 3:** The images in a voucher flora must have the current taxon name to the top right of the  
652 species images (Figure 11). Record taxonomic concept references, name changes in client  
653 voucher flora, and additional notes below the image (Figure 11). Specifically, the taxonomic  
654 concept should be referenced with details (e.g., plate, figures) from publications (see example  
655 Figure 12C,E), or the plate name and OTU from publicly available voucher flora (see example  
656 Figure 12A,D).

657

Taxon name



\*Taxonomic concept reference

\*\*Notes regarding name changes in client voucher

\*\*\*Additional helpful notes

658

659 **Figure 11.** *Labeling and notes template for a taxon in a voucher flora plate.*

660 **Step 4:** Email all phycologists when the project voucher flora is complete, and request reviews  
661 by at least two.

662 **Step 5:** Review comments and requests for additional raw images to aid in a further review  
663 should be sent via email. A meeting to discuss opinions is held to evaluate all the review  
664 comments.

665 **Step 6:** Once a consensus is reached for a new taxon it can be added to the on-going client  
666 voucher flora. The voucher flora is consulted and cited (e.g., NEON\_Achnanthidium\_plate1.psd)  
667 for future identifications. Consultation with external phycologists may be required if a  
668 consensus cannot be reached on the naming of a new taxon.

669 **Step 7:** Create a “species complex” if necessary. A complex or “slash group” may be warranted  
670 if phycologists cannot agree on how to separate morphotypes consistently (e.g., *Achnantheidium*  
671 *minutissimum* complex, Figure 12A).

672 **Step 8:** Export the finished project voucher flora in two formats: psd and jpg. A GIMP specific  
673 format “.psd” will allow users to copy/paste images into the on-going client voucher flora. The  
674 shareable static version “.jpg” will allow voucher flora plates to be delivered easily to clients  
675 upon completion of the project.

676


### 677 ***Updates to client voucher flora***

678 Updates to species concepts in the client voucher flora are often needed, given the ever-  
679 evolving field of diatom taxonomy and interpretation of species concepts. It is essential that the  
680 timing/date of concept changes be documented in the on-going client vouchers. Several  
681 situations that lead to revision of concepts in client voucher floras include:

- 682 • Name change in the literature (i.e., a name change based on the same type specimen, see S.  
683 Spaulding DONA post on 23 May 2023)(Spaulding 2023a; b). Cited literature can be a recent  
684 publication or new literature discovered by Rhithron phycologists. A descriptive note should  
685 be made in the client voucher beneath the taxon images (e.g., 2010-2017 reported as  
686 *Diatoma mesodon*, updated 27 February 2018 to *Odontidium mesodon*, Figure 12B).
- 687 • A newly described taxon is published which clearly resembles an existing taxon in the client  
688 voucher flora. A descriptive note should be made in the client voucher beneath the images  
689 (e.g., 2014-2020 reported as *Gomphosphenia* sp 1, updated 12 January 2021 to  
690 *Gomphosphenia praegnans*, Figure 12E).
- 691 • Two species (or more) in the client voucher are deemed “problematic” for consistent  
692 separation and need to be combined. Typically, the problematic species do not have solid  
693 distinguishing features visible with DIC optics that are resolvable among the phycologists. The  
694 issue often becomes apparent during review of internal QCs and subsequent years of  
695 reviewing and using the client voucher. A descriptive note should be made in the client  
696 voucher beneath the combined images (e.g., 2008-2010 reported separately as *G. parvulum*  
697 and *G. exillissimum*, combined 6 May 2011 into *Gomphonema parvulum* complex, Figure  
698 12D).
- 699 • A single species needs to be split into 2 species (or more) in the client voucher due to the  
700 realization that solid distinguishing features are apparent and agreed upon by the  
701 phycologists. Often this realization comes about due to a recently published or discovered  
702 reference that offers clarity and guidance. A descriptive note should be made in the client  
703 voucher beneath the images of the taxa that were split out (e.g., 2009-2017 reported as  
704 *Achnantheidium minutissimum*, updated 3 March 2018 to *Achnantheidium affine*, Figure 12C).

705

706



707

708

709

**Figure 12.** Example of an on-going client voucher with taxonomic concept references and the types of descriptive notes for recording name/nomenclature changes.

710

711

712

### 713 **Metadata for image reference collections**

714 Additional metadata is documented for all new taxa identified in samples. Metadata for a new  
715 taxon includes image location, image magnification, XY coordinates of a representative  
716 specimen, and taxonomic references.

717

718 **Step 1:** Add metadata for a new taxon entry by clicking on the plus symbol under the Reference  
719 column in EPIC. All metadata recorded in steps 2-5 need to be entered into this Reference  
720 Collection tab (Figure 13).

721 **Step 2:** Record the magnification used to take the diatom image (typically 1000X).

722 **Step 3:** In the Location (Vernier coordinates) field, record the Vernier (XY) coordinates to  
723 pinpoint the specimen location on the slide. **Record XY coordinates immediately when**  
724 **encountering a new taxon for the client voucher** as it is time-consuming to “back track” to the  
725 specimen location.

726 **Step 4:** Record the taxonomic reference used for identification in the Comment field (Figure  
727 13). Include the full citation information available for the reference (e.g., Potapova, M. (2009).  
728 Achnantheidium minutissimum. In Diatoms of North America. Retrieved February 19, 2024, from

729 [https://diatoms.org/species/achnanthidium\\_minutissimum](https://diatoms.org/species/achnanthidium_minutissimum)). Commonly referenced USGS  
730 voucher floras listed on DONA's voucher flora archive webpage  
731 (<https://diatoms.org/practitioners/voucher-flora-archive>) (Spaulding *et al.* 2021) can be cited  
732 with abbreviations (SESQA, NESQA, PNSQA), however include plate number (e.g.,  
733 gomphonema1\_MERGEreconcile\_PNSQA) and the OTU code (e.g., PNSQA GOM52).

734  
735 **Step 5:** In the Filename field, record the location of the reference image. Images in voucher  
736 floras will have the name of the specific project plate where the image/taxon is located (e.g.,  
737 NEON24\_Achnanthidium\_plate1.jpg).

738

The screenshot shows a dialog box titled "Add to Reference Collection" with a close button (X) in the top right corner. The species name "Achnanthidium minutissimum" is displayed at the top. Below it are several input fields: "Filename" with the value "NEON24\_Achnanthidium\_p", "Magnification:" with the value "1000", and "Location (Vernier coordinates):" with "X" set to "38" and "Y" set to "107". A "Comments:" section contains a text area with the following text: "Potapova, M. (2009). Achnanthidium minutissimum. In Diatoms of North America. Retrieved February 16, 2024, from https://diatoms.org/species/achnanthidium\_minutissimum". At the bottom of the dialog are "Save" and "Close" buttons.

739  
740 **Figure 13.** Example of metadata entered for a new taxon.  
741

742

## 743 Quality Assurance and Quality Control

### 744 QA/QC summary

745 To ensure project integrity, these important rules are followed:

- 746 • Phycologists work as a team and continually communicate to improve data quality for the  
747 entire project.
- 748 • Phycologists adhere to the client-specified protocol (e.g., target valve count, taxonomic  
749 resolution, new taxa/image documentation).
- 750 • Project images and voucher flora are reviewed by all Rhithron phycologists, to assure that  
751 identifications are accurate and consistent.

- 752 • All client-requested deliverables are provided, including citations for reference materials and  
753 project/client voucher floras.
- 754 • All calculations are double-checked for accuracy.
- 755 • A summary of QA/QC procedures and results are documented and delivered to the client  
756 along with a summary of sample processing procedures.
- 757 • Re-identification and enumeration of the same transect(s) is conducted for at least 10% of  
758 samples to determine comparability of sample analysis, such as the Percent Taxonomic  
759 Disagreement parameter (PTD), (Stribling *et al.* 2003, 2008)). Rhithron follows the NRSA 2023-  
760 2024 QAPP minimum data quality objective (MDQO), which is  $PTD \leq 25\%$  (USEPA 2022).

761

### 762 ***Internal re-identification of identified transect(s)***

763 Ten percent of samples from a client’s project are identified by a second Rhithron phycologist  
764 along the same scribed transect line(s) that were used in the original identification. QC samples  
765 are randomly selected, and the re-identification proceeds as a “blind” examination, that is, the  
766 QC phycologist does not have access to the original identification data. The initial identifications  
767 are compared with the QC identifications, and the Percent Taxonomic Disagreement (Stribling  
768 *et al.* 2003) parameter is calculated (Figure 14).

769 Percent Taxonomic Disagreement is calculated as:

$$770 \quad PTD = \left(1 - \left(\frac{comp_{pos}}{N}\right)\right) \times 100$$

771 Where  $comp_{pos}$  is the number of agreements and  $N$  is the total number of taxa identified.

772

773 Additionally, the relative percent difference in transect lengths (RPDTL, Figure 14) is calculated  
774 to determine if the same transects were read. No MDQO exists for the relative percent  
775 difference in transect length parameter. Taxonomic and transect length discrepancies are  
776 examined and discussed between the original phycologist and the QC phycologist.

777 Relative Percent Difference of Transect Length (RPDTL) is calculated as:

$$778 \quad RPDTL = \left(\frac{(TL_1 - TL_2)}{\frac{(TL_1 + TL_2)}{2}}\right) \times 100$$

779 Where  $TL_1$  and  $TL_2$  are the Transect Lengths in mm reported from the primary and secondary  
780 taxonomists respectively.

$$781 \quad PDE = \left(\frac{Count_1 - Count_2}{Count_1 + Count_2}\right) \times 100$$

782 Where Count 1 and Count 2 are the Sum total of all diatoms enumerated by each phycologist.

783

784

785

786

$$PTC_{abs} = ABS \left( \left( \frac{Target_{count1}}{Total_{count1}} \right) - \left( \frac{Target_{count2}}{Total_{count2}} \right) \right)$$

787 Where the Target counts are the number of individuals identified to the required taxonomic  
788 target (usually species) for each respective taxonomist, and Total Counts are the total number  
789 of individuals identified in the sample from each respective taxonomist.

790

791 Difficult taxa are further discussed with the entire phycology team to come to a consensus on  
792 the handling of the taxa, and to provide a basis for taxonomic precision. Data changes or  
793 nomenclatural edits are applied to the entire project dataset, not just the sample that  
794 underwent QC. Changes to taxonomic concepts in the client voucher flora may also be applied  
795 warranted.

796

BA23D01001			Diatoms					
			Analysis 1	Analysis 2	Quality Control Statistics			
			Original		Rectified			
TOTAL COUNT:			601	600	PDE	0.08%	PDE	0.08%
Identification Date:			12/19/2023	12/11/2023	RPDE	0.04%	RPDE	0.04%
Taxonomist:			JT	LV	PTD	9%	PTD	9%
TRANSECT LENGTH			8.2	7.5	B-C	91%	R(B-C)	91%
RPD Transect Length (%)			4.46%					
			Original Counts			Rectified Counts		
Taxon	Original Comments	QC Comments	W	A	B			QC Rectifying Comments
Achnanthisidium	Achnanthisidium sp8 SESQA		2	0	2	0	2	0
Achnanthisidium minutissimum			460	477	460	460	477	460
Amphora pediculus			0	1	0	0	1	0
Craticula accommoda			2	1	1	2	1	2
Craticula molestiformis			0	4	0	0	4	0
Craticula sp. 3 SESQA			5	9	5	5	9	5
Encyonema minutum			13	15	13	13	15	13
Eunotia	girdle	girdle view	2	2	2	2	2	2
Eunotia minor			1	1	1	1	1	1
Fragilaria	small form, Fragilaria sp.23 SESQA?		2	0	2	0	2	0
Fragilaria voucheriae			0	2	0	0	2	0
Gomphonema auritum		3 girdle view added	3	12	3	12	3	12
Gomphonema lagenula		7 girdle view added	39	27	27	39	27	39
Gomphonema parvulum			8	3	3	8	3	8
Mayamaea atomus			4	2	2	4	2	4
Mayamaea cahabaensis			1	0	0	1	0	1
Navicula cryptocephala			4	2	2	4	2	4
Navicula gregaria			0	1	0	0	1	0
Navicula microcari			0	1	0	0	1	0
Nitzschia elegantula	mostly girdle		6	0	0	6	0	6
Nitzschia frustulum	narrow form SESQA	7 girdle view added	38	27	27	38	27	38
Nitzschia palea			0	5	0	0	5	0
Nitzschia palea v. debilis			8	2	2	8	2	8
Nitzschia sp. 34 SESQA			0	2	0	0	2	0
Pinnularia	obscured		1	0	0	1	0	1
Sellaphora	obscured, S.nigri?		1	1	1	1	1	1
Sellaphora pupula	small form in voucher		1	0	0	1	0	1
Staurisira construens v. venter			0	3	0	0	3	0

797

798 **Figure 14.** Internal diatom QC analysis example for re-identification of same transect(s).  
799 Calculations are made in Excel within locked cells in order to avoid miscalculation of QC metrics.  
800 Important fields to note in the spreadsheet include the Analysis 1 and Analysis 2 fields, and the  
801 Quality Control Statistics fields. Analysis 1 includes the counts made in the original  
802 identification, and Analysis 2 includes the counts made by the QC phycologist. Where “0” occurs  
803 in one of these fields, a discrepancy in the identification of one or more specimens is indicated.  
804 In the Quality Control Statistics fields, calculations for the Percent Difference in Enumeration  
805 (PDE), and the Percent Taxonomic Disagreement (PTD), among other parameters, are given.  
806 Values for these parameters are given for a preliminary comparison of the two analyses in the  
807 Original (green filled) fields. During the QC process, taxonomic discrepancies are discussed.  
808 Discrepancies may be attributable to changes in nomenclature, further identification of difficult  
809 specimens, or duplication of entries, etc. If consensus is reached, the data are rectified to reflect  
810 the agreed-upon changes, which are elaborated in the Taxon and Comments fields. A

811 *comparison of the rectified data is then conducted, and quality control statistics are computed*  
 812 *and reported in the Rectified (red filled) fields.*

813 **Table 3: Quality control parameters, associated measurement quality objectives, and potential corrective**  
 814 **actions.**

Quality Control Parameter	Measurement Quality Objective	Corrective Action
Percent Taxonomic Disagreement (PTD)	25%	Review of taxa of concern, make necessary data changes throughout sample set. Additional samples may need secondary identification. Voucher adjustments
Percent Difference in Enumeration (PDE)	1%	Discuss slides, optics and counting requirements review SOP- to assure that counting technique is faithfully followed.
Relative Percent Difference in Transect Length (RPDTL)	10%	Discuss slides, optics, and counting, External QC may result in higher RPDTL due to FOV widths among scopes
Percent Taxonomic Completeness Absolute Difference (PTCabs)	10%	Review of voucher concepts, discuss obscured views, or girdle views

815

816 ***External re-identification of identified transect(s) by external phycologist***  
 817 **NEON WORKING GROUP NEEDS TO MAKE A DECISION ON WHETHER TO INCLUDE RE-**  
 818 **ANALYSIS OF XX% SLIDES IN PROJECTS BY AN EXTERNAL/NON-RHITHRON PHYCOLOGIST.**  
 819 **ADDITIONALLY, OR ALTERNATIVELY, AN EXTERNAL REVIEW OF PROJECT/CLIENT VOUCHERS**  
 820 **COULD BE DONE PERIODICALLY.**

821 Rhithron will send analyzed diatom slides to another phycologist for external QC upon client  
 822 request. Pertinent information will be shared with the external phycologist so that, ideally, the  
 823 re-analysis can be conducted exactly like Rhithron’s internal QC procedures (e.g., same  
 824 transects etc.). The project and on-going client voucher flora will be provided to the external  
 825 QC taxonomist to review and use during re-analysis of the QC slides. Experience has shown that  
 826 common issues that arise with external QC comparisons are 1) original transects are not  
 827 located/followed and 2) microscopes have different field of view dimensions (note: all Rhithron  
 828 phycologist scopes have 250-micron diameter FOVs at 1000X). These two issues complicate  
 829 diatom QC comparisons and are documented when comparing the taxa lists and calculated  
 830 similarity metrics. Rhithron phycologists review and discuss the results as a group. A summary

831 presentation is drafted during the discussion that highlights QC metric scores, images of  
832 difficult/problematic taxa and potential actions to rectify the project data. Subsequently, the  
833 external QC summary is shared with the external phycologist. The summary is discussed to  
834 guide any changes that need to be made to project data and vouchers.


835

## 836 **Data Review**

837 Diatom data for each project undergoes a multistep review procedure, which results in  
838 thorough reviews and data corrections (Figure 15). Initial and second reviews are completed,  
839 and the data is sent on to Rhithron’s Quality Assurance (QA) Officer. In addition to data review,  
840 the QA Officer highlights any new taxa in the project that were not previously on the client’s  
841 taxa list or included in the voucher. The initial phycologist makes corrections and edits the data  
842 after each review. Common edits to diatom datasets include:

- 843 • Genus-level taxon entries that need a description in Comment field regarding why target-level  
844 was not reported. Description examples might include girdle views only, obscured valves,  
845 damaged, and provisional names (e.g., *Nitzschia* sp1).
- 846 • Nomenclature needs to be brought up to date.
- 847 • Metadata is missing or does not agree with the Project Inventory Report benchsheet. Project  
848 metadata, especially those data transcribed from benchseets to data fields, are evaluated for  
849 accuracy, precision and completeness during the initial data review by the original  
850 phycologist, secondary phycologists and the QC officer.
- 851 • Duplicate entries.
- 852 • Total count less than target without comment “maximum effort 5 transects/4h”.
- 853 • Transect length discrepancy in Sample Comments.
- 854 • Review bench sheet entries that match data entry fields.

855



856  
857  
858

**Figure 15.** *Data review flow chart.*

## 859 **Instrument and equipment inspection, maintenance and calibration**

860 All microscopes and equipment needed to process diatom samples are inspected regularly  
861 according to manufacture specifications. Rhithron contracts with Pacific Microsystems LLC,  
862 Seattle, Washington to perform annual maintenance on all microscopes. Digital image  
863 calibrations are routinely checked against stage micrometers to assess for potential drift in  
864 measurement accuracy.

865

## 866 **Data Delivery and Sample Material Shipment**

867 Finalized data are transferred from EPIC to RAILIS by the Data Specialist via an electronic  
868 upload. All data, including taxonomic data and metadata are formatted to the client's  
869 specifications and uploaded to the NEON data repository. Labeled and barcoded material are

870 retrieved from the phycologists, sample material barcodes are scanned and additional  
871 metadata from RAILIS are formatted to the necessary shipping manifest specifications. Sample  
872 shipping manifests are uploaded to NEON and sample materials are securely packaged and sent  
873 to the biorepository (ASU). Sample fractions will be in three fractions 1) air-dried 20mL glass  
874 scintillation vial of digested diatom material, 2) slide that was used for the identification and  
875 enumeration, and 3) the duplicate slide that is not associated with any taxa count data.  
876 Scintillation vials will be packaged in flats of 100 (10x10) with cardboard dividers. Slides will be  
877 packaged in 100 count slide boxes (50 samples per box), additional padding will be added to  
878 slide boxes to assure minimal movement of slides during transport.

879

## 880 **Accrual Year Data Review**

881 After all samples from NEON's designated accrual year have been successfully processed a  
882 cumulative data review will commence. While similar to the project data reviews, special  
883 attention is paid to the voucher flora and consistency over the accrual year's samples. During  
884 this data review taxa are discussed from throughout the accrual year, and any previous years,  
885 to assure precision among sample sets. Any necessary data changes to sample taxa count  
886 information or vouchers are communicated to NEON staff. The image voucher is finalized and  
887 an associated taxa translation table (alias tracking) spreadsheet containing synoptic notes is  
888 created to accompany the version-controlled Voucher Flora. The new finalized voucher flora in  
889 both PDF and .psd formats are delivered to NEON biologists for archive.

890

891  
892  
893  
894  
895  
896  
897  
898  
899  
900  
901  
902  
903  
904  
905  
906  
907  
908  
909  
910  
911  
912  
913  
914  
915  
916  
917  
918  
919  
920  
921  
922  
923  
924  
925  
926  
927  
928  
929  
930  
931

## References

Bishop, I.W., Espisito, R.M., Tyree, M. & Spaulding, S.A. (2017) A diatom voucher flora from selected southeast rivers (USA). *Phytotaxa* 332, 101.  
<https://doi.org/10.11646/phytotaxa.332.2.1>

Charles, D.F., Knowles, C. & Davis, R.S. (2002) Protocols for the analysis of algal samples collected as part of the U.S. Geological Survey National Water-Quality Assessment Program. *Report No. 02-06*, 124.

Lee, S.S., Spaulding, S.A., Stevenson, R.J., Edlund, M.B., Heinlein, J. & Burge, D. (2020) Webinar:Voucher flora workshop. *Diatoms of North America*. Available from:  
<https://diatoms.org/news/online-voucher-flora-workshop>

Potapova, M.G., Lee, S.S., Spaulding, S.A. & Schulte, N.O. (2022) A harmonized dataset of sediment diatoms from hundreds of lakes in the northeastern United States. *Scientific Data* 9, 1–9. <https://doi.org/10.1038/s41597-022-01661-3>

Rhithron (2014) EPIC Instruction Manual for Rhithron Phycologists. , 10.

Rhithron (2021) Periphyton Standard Operating Procedure ver.21.2.b. , 15.

Rhithron (2022) Health and Safety Manual ver. 2.0. , 11.

Spaulding, S.A. (2023a) FAQ: What are homotypic and heterotypic synonyms? Available from:  
<https://diatoms.org/news/faq-what-are-homotypic-and-heterotypic-synonyms>

Spaulding, S.A. (2023b) FAQ: What is a basionym? *Diatoms of North America*. Available from:  
<https://diatoms.org/news/faq-what-is-a-basionym>

Spaulding, S.A., Potapova, M.G., Bishop, I.W., Lee, S.S., Gasperak, T.S., Jovanoska, E., Furey, P.C. & Edlund, M.B. (2021) Diatoms.org: supporting taxonomists, connecting communities. *Diatom Research*. <https://doi.org/10.1080/0269249X.2021.2006790>

Stribling, J.B., Moulton II, S.R. & Lester, G.T. (2003) Determining the quality of taxonomic data. *Journal of the North American Benthological Society* 22, 621–631.  
<https://doi.org/10.2307/1468357>

Stribling, J.B., Pavlik, K.L., Holdsworth, S.M. & Leppo, E.W. (2008) Data quality, performance, and uncertainty in taxonomic identification for biological assessments. *Journal of the North American Benthological Society* 27, 906–919. <https://doi.org/10.1899/07-175.1>

Tyree, M.A., Carlisle, D.M. & Spaulding, S.A. (2020) Diatom enumeration method influences biological assessments of southeastern usa streams. *Freshwater Science* 39.  
<https://doi.org/10.1086/707725>

Tyree, M.A., Garrison, P., Enache, M. & Spaulding, S.A. (2018) A diatom voucher flora of lakes from the Northeast United States:Voucher Edition. , 155.

USEPA (2017) National Rivers and Streams Assessment 2018-19: Laboratory Operations Manual.EPA-841-B-17-004. , 185.

USEPA (2022) National Rivers and Streams Assessment 2023/24: Quality Assurance Project Plan.EPA-841-B-22-004. , 119.