

Standard Operating Procedure: Receiving NEON KCl Extracts

5/18/2017

Prepared by D. Reuss, Lead Analyst
EcoCore Laboratory, Colorado State University

1. Take the cooler to the bench near the NEON freezer.
2. Open the cooler and step away for a minute to allow the carbon dioxide (CO₂) to dissipate. CO₂ is an asphyxiation hazard so move to fresh air and get help if you feel dizzy or light headed.
3. Remove the sample manifest.
4. Remove enough of the packaging material and dry ice to be able to access the samples. Wear gloves when handling dry ice!
5. Note the temperature state of the samples and write it on the manifest along with the date and your name. The temperature states are: Frozen, Thawed but cold, or Warm.
6. Check the sample labels on the vials against the manifest. Make a note on the manifest of any missing samples, illegible or missing labels, broken or damaged vials or caps, and any signs of leaks.
7. Put the samples in the freezer promptly after inspection.
8. Add the samples to the monthly task log including the temperature state and any problems with the sample condition.
9. Notify NEON if the samples were warm or there were any problems with the sample condition.
10. Put the manifest in the NEON file drawer in room A250.

Standard Operating Procedure: Preparing Reagents for NO₃, and NH₄ analysis

5/26/2017

Prepared by D. Reuss, Lead Analyst
EcoCore Laboratory, Colorado State University

Calibration standards and check standards should both be prepared from primary standard grade salts or other certified or traceable materials, but the two should come from entirely different sources to be able to use them to verify each other.

Materials and Cleaning

All solid standard materials (salts) are oven dried at 105° C overnight, then placed in desiccators to cool to room temperature before weighing. All weighing is performed on a 4 place (0.0001 g) balance using a clean spatula and a disposable plastic weighing boat. Balances are calibrated by a certified balance technician annually and performance is verified using check weights periodically by the lab managers.

All solutions are made in Class A glass volumetric flasks and all dilutions are made using Class A glass volumetric pipettes and volumetric flasks. Standards are stored in Nalgene polyethylene or polypropylene bottles. All calibration and check standards for NEON KCl ammonium analyses are prepared using 2 molar KCl. The KCl solution is prepared using DI water with an electrical conductivity of less than 1 micro Siemen. All calibration and check standards for NEON KCl nitrate analyses are prepared using ultrapure DI water as above.

All storage bottles, volumetric flasks and pipettes, funnels and beakers used for standard preparation are acid washed prior to use. The acid washing procedure is: items are individually rinsed 3 times with DI water. Items are then fully immersed in a 1.2 Molar Hydrochloric Acid bath (1 part concentrated HCl to 9 parts DI water) and left overnight. Items are then individually rinsed inside and out 6 times with DI water. Volumetric flasks are stored full of DI water and emptied just before use. Volumetric pipettes are oven dried at 105° C. Beakers and funnels and storage bottles are air dried. All materials are stored in a designated space separate from other lab glassware.

Calibration and check standards are prepared using primary standard grade salts. Ammonium standards are prepared using Ammonium Sulfate, nitrate standards are prepared using Potassium Nitrate and the nitrite check standard is prepared using Potassium Nitrite.

Calibration Standards

A series of calibration standards are prepared and analyzed to create a standard curve which is used to calibrate the instrument. The concentrations of these standards are chosen to cover the expected range of the samples in milligrams nitrogen per liter. 5-9 standards are typically run in a calibration curve. A significant volume (typically 500 ml) of a stock solution with a high concentration (e.g. 1000 mg-N/l) is prepared. This allows for the amount of salt to be weighed to fall in the 2-5 gram range for both analytes. The salt is weighed to +/- 0.0001 grams of the nominal weight on a 4 place analytical balance. This allows for an extremely small weighing error. Stock solutions are prepared in DI water.

Working standards are prepared by diluting the stock solution with 2 M KCl (for ammonium) and DI water (for nitrate) using volumetric pipettes and flasks. Very low concentration standards (less than 0.5

ppm) are prepared by a two-step serial dilution. The stock solution is diluted to 10 or 1 ppm, then this solution is diluted to make the final working standards.

Check Standards

Check standards are prepared using the same techniques specified above for the working standards, but using a separate chemical source.

Storage

All calibration standards, check standards and stock solutions are stored in a refrigerator from day to day and frozen if storage is for more than a week. Stock solutions may be stored up to one year if frozen. Working standards may be stored up to 6 months refrigerated.

Ammonium Reagents

NH₄ Buffer

Dissolve 24 grams sodium hydroxide, 58 grams sodium citrate, and 50 grams disodium EDTA in approximately 800 ml DI water in a 1 liter volumetric flask. Add DI water to the 1 liter volume. Filter through a 10 micron filter. May be stored at room temperature for up to 4 months. Add 6 drops Brij -36 solution to 500 ml buffer before use on the instrument. This working buffer is good for up to 2 weeks.

Salicylate- Nitroferricyanide

Dissolve 300 grams Sodium Salicylate and 0.6 grams sodium nitroferricyanide in approximately 400 ml DI water in a 1 liter volumetric flask. Add DI water to the 1 liter volume. Filter through a 10 micron filter. Store refrigerated for up to 2 months.

Sodium Hypochlorite

Add 12 ml household bleach (5.25% sodium Hypochlorite) to 188 ml. DI water. Make fresh daily.

Nitrate Reagents

Ammonium Chloride Buffer

Dissolve 638 grams of ammonium chloride and 0.75 grams of disodium EDTA in approximately 6 liters of DI water. Add DI water to 7.5 liters total volume. Adjust to pH 8.5 +/- 0.1 using ammonium hydroxide. Store at room temperature for up to 3 months.

NO₃ Color Reagent

Slowly add 100 ml concentrated phosphoric acid (Caution: severe contact hazard) to approximately 800 ml DI water in a 1 liter volumetric flask. Dissolve 40 grams of sulfanilamide and 2 grams N-1-naphthylethylenediamine in the diluted phosphoric acid and add DI water to the 1 liter volume. Filter through a 10 micron filter. Store refrigerated for up to 2 months.

2% Cupric Sulfate

Dissolve 2 grams cupric sulfate in 100 ml DI water. Store at room temperature for up to 5 years.

0.2% Cupric Sulfate

Add 10 ml of 2% Cupric Sulfate to 90mls. DI water. Store at room temp

Standard Operating Procedure: Setting up and running NH₄/NO₃ analyses on the Alpkem

5/21/2018

Prepared by Guy Beresford
EcoCore Laboratory, Colorado State University

Initial Set-up

1. Remove NO₃ color reagent and NH₄ salicylate-nitroferricyanide bottles from the fridge. Place the NO₃ color reagent in a warm water bath to warm up to room temperature.
2. Check sampler wash waste reservoir under bench to right of instrument. Empty into sink if more than half full.
3. Turn on main power switches for pump, main unit and wash reservoir.
4. Open FlowView software.
5. Lock down the pump platens and pull the engaging levers straight up for all the pump tubes labeled NH₄ and NO₃ plus “to sampler wash” and “from sampler wash”.
6. Press the “run stop” button on the pump, the display should read -50.0.
7. Empty, rinse and refill the large DI water bottle for the sampler wash and connect line.
8. Empty, rinse and refill a large DI water bottle for the NO₃ carrier, color and buffer lines and connect.
9. Empty, rinse and refill a second large DI water bottle
10. Connect the NH₄ lines: buffer, hypochlorite, and salicylate and carrier lines to the second DI wash bottle.
11. Let pump for 10 minutes and check for leaks.
12. Discard old hypochlorite solution and make fresh. 12 mls bleach and 178 mls DI water.
13. Move the NO₃ buffer line to the buffer bottle.

Condition the cadmium column

1. Place the column, syringe with tube to fit column, open specimen cup of NO₃ buffer, and open cup of 0.2% CuSO₄ on the bench within easy reach.
2. Draw about 1 ml of buffer into the syringe and remove all air bubbles from the tube.

3. Break the tubing connection on the column and elevate the connector end slightly until a drop of liquid appears at the tube end.
4. Hold the syringe plunger end up and depress the plunger until a drop of buffer appears at the connector end.
5. Connect the syringe to the column without introducing any air.
6. Place the free end of the column into the buffer cup and slowly depress the plunger to expel the air from the free end of the column.
7. Transfer the free end to the cup of CuSO_4 and draw up about 0.5 ml.
8. Depress the plunger slightly to release any tension on the liquid in the column and transfer the end to the cup of buffer.
9. Draw up 3-5 mls of buffer, and then depress the plunger slightly to release any tension on the liquid in the column.
10. Remove the end from the buffer, place the column, tubing and syringe flat on the bench, then disconnect the syringe and connect the two ends of the column together.
11. Take the column and syringe to the sink and rinse the outside with DI water. Put the buffer and CuSO_4 cup away and wipe the bench with a wet sponge. The NO_3 buffer contains NH_4 , which can contaminate your samples!

Connect the cadmium column to the instrument

1. Make sure there are no bubbles in the nitrate cartridge (about 10 minutes after moving the NO_3 buffer line to the buffer bottle).
2. Press the "run stop" button on the pump, press the "mode" button 4 times, the display should read "2.0", press the "run stop" button again, this puts the pump speed to minimum.
3. Break the connection in the line between the "to column" and "from column" ports on the NO_3 cartridge.
4. Break the tubing connection on the column and note which end will connect to the line from the "to column" port on the cartridge. Elevate the other end of the tubing until a drop of liquid forms at the end. Wait until a drop forms on the end of the "to column" line and connect it to the column without introducing any air.
5. Connect the other end of the column to the "from column" line on the cartridge.
6. Press the "run stop" button on the pump, press the "mode" button 3 times, the display should read "-50.0", press the "run stop" button again, this puts the pump speed back to normal.
7. Move the NO_3 color reagent line to the color reagent bottle.

8. Move the NH₄ reagent lines to their respective reagent bottles.
9. Make sure the line from the sampler is completely clear of all other lines and the reagent bottles so the sample probe can move freely.
10. Let the system pump at least 15 minutes before beginning a sample run.

Set up the sample table

1. Press the sample table button on the main tool bar in FlowView. Alternatively, import samples names into blank sample table Excel csv file.
2. Type in "cup #", "name" and "type" for each sample and standard to be analyzed. Refer to the typical sample table example in *SOP: Setting up the sample table for NEON KCl* (pg 9).

Starting a sample run

1. Click the configure analyzer icon and then the system icon. Enable chemistries to used and click save.
2. Click home icon and Run Samples icon. To right of sample table browse for sample table file and select. Enter operator ID and result file name. Click on start.
3. The data collection window will appear on the screen. At this point, the software will monitor the baselines for 60 minutes and then start sampling.
4. Debubble both flow cells by pinching and releasing the outflow tubing. Repeat until no bubbles appear.
5. Pour standards and samples into the cups and place in the sampler. Get at least 20 cups poured before starting the sampler.
6. Check baselines for stability and drift, should be less than 500 micro-absorbance units.
7. Press the Start icon to begin sampling.
8. The first NO₃ peak will appear about 1 minute after the first sample is drawn. The first NH₄ peak will appear in about 4 minutes.

Instrument shutdown

1. Press the "run stop" button on the pump to stop the reagent flow.
2. Disconnect the cadmium column from the NO₃ cartridge and connect the two ends together. Be careful to minimize the amount of air let into the column.

3. Connect the “to column” and “from column” lines together.
4. Press the “run stop” button on the pump to restart the flow.
5. Move the NO₃ lines to the NO₃ rinse bottle and the NH₄ lines to the NH₄ manifold startup solution bottle. Let pump about 10 minutes.
6. Disconnect all the lines and let the system pump air for about 15 minutes or until all the liquid is out of the instrument.
7. Press the “run stop” button to stop the pump and turn off the main power switch to the instrument.
8. Push the engaging levers on all the pump platens to horizontal and disconnect one end of each platen from the pump.
9. Empty the NH₄ waste into the container on the floor next to the instrument.
10. Pour the NO₃ waste into the container on the floor next to the instrument.
11. Empty and discard the sample cups and tubes. Rinse the sampler racks with DI water and set on the bench to dry.
12. Return the NO₃ color reagent and NH₄ salicylate-nitroferricyanide bottles to the fridge.
13. Refrigerate the standards if the instrument will be used the next day or freeze them if it will be longer.
14. Wipe down the instrument and benches with a moist sponge.

Standard Operating Procedure: Setting up the sample table for NEON KCl

10/08/2018

Prepared by G. Beresford

EcoCore Analytical Services, Colorado State University

This protocol covers Sample Table setup to ensure proper QA/QC for NEON KCl extracts. KCl extracts are run twice. Standards and checks are made in DI water for the NO₃ run, with the reagent blank also being water. For the NH₄ run, standards and checks are made in KCl and the reagent blank is also KCl.

Creating the Sample Table

1. Set up the Alpkem for NH₄/NO₃ analysis according to the instrument directions above.
2. Copy NEON sample IDs into an Excel .csv sample table file, following these examples.

Nitrate sample table example

908	Sync DI	1 SPL	Cup1	1						1	0,
909	10 NO2 ch	1 SPL	Cup2	1						1	0,
	0 blank	1 SPL	Cup2	1						1	0,
	0 blank	1 SPL	Cup2	1						1	0,
	0 blank	1 SPL	Cup2	1						1	0,
901	baseline	1 RB	Cup2	1						1	0,
901	0	1 STD1	Cup1	1						1	0,
902	0.05	1 STD2	Cup1	1						1	0,
903	0.1	1 STD3	Cup1	1						1	0,
904	0.5	1 STD4	Cup1	1						1	0,
905	1	1 STD5	Cup1	1						1	0,
906	2	1 STD6	Cup1	1						1	0,
907	5	1 STD7	Cup1	1						1	0,
908	10	1 STD8	Cup1	1						1	0,
	0 blank	1 SPL	Cup2	1						1	0,
	0 blank	1 SPL	Cup2	1						1	0,
901	baseline	1 RB	Cup1	1						1	0,
910	10 10 chec	1 ICV	Cup1	1						1	0,
	0 blank	1 SPL	Cup2	1						1	0,
901	baseline	1 RB	Cup1	1						1	0,
	0 blank	1 SPL	Cup2	1						1	0,
909	10 NO2 ch	1 SPL	Cup1	1						1	0,
	0 blank	1 SPL	Cup2	1						1	0,
901	baseline	1 RB	Cup2	1						1	0,
101		1 SPL	Cup2	1						1	0,
102		1 SPL	Cup2	1						1	0,
103		1 SPL	Cup2	1						1	0,
104		1 SPL	Cup2	1						1	0,
105		1 SPL	Cup2	1						1	0,
106		1 SPL	Cup2	1						1	0,
107		1 SPL	Cup2	1						1	0,
108		1 SPL	Cup2	1						1	0,
109		1 SPL	Cup2	1						1	0,
110		1 SPL	Cup2	1						1	0,
901	baseline	1 RB	Cup1	1						1	0,
910	10 10 chec	1 ICV	Cup1	1						1	0,
	0 blank	1 SPL	Cup2	1						1	0,
901	baseline	1 RB	Cup1	1						1	0,
111		1 SPL	Cup2	1						1	0,
112		1 SPL	Cup2	1						1	0,

Ammonium sample table example

908	Sync KCl	1 SPL	Cup1	1						1	0,
909	10 NO2 ch	1 SPL	Cup1	1						1	0,
	0 blank	1 SPL	Cup2	1						1	0,
	0 blank	1 SPL	Cup2	1						1	0,
901	baseline	1 RB	Cup2	1						1	0,
901	0	1 STD1	Cup1	1						1	0,
903	0.1	1 STD3	Cup1	1						1	0,
904	0.5	1 STD4	Cup1	1						1	0,
905	1	1 STD5	Cup1	1						1	0,
906	2	1 STD6	Cup1	1						1	0,
907	5	1 STD7	Cup1	1						1	0,
908	10	1 STD8	Cup1	1						1	0,
	0 blank	1 SPL	Cup2	1						1	0,
909	10 NO2 ch	1 SPL	Cup1	1						1	0,
	0 blank	1 SPL	Cup2	1						1	0,
901	baseline	1 RB	Cup1	1						1	0,
910	5 5 check	1 ICV	Cup1	1						1	0,
	0 blank	1 SPL	Cup2	1						1	0,
901	baseline	1 RB	Cup1	1						1	0,
101		1 SPL	Cup2	1						1	0,
102		1 SPL	Cup2	1						1	0,
103		1 SPL	Cup2	1						1	0,
104		1 SPL	Cup2	1						1	0,
105		1 SPL	Cup2	1						1	0,
106		1 SPL	Cup2	1						1	0,
107		1 SPL	Cup2	1						1	0,
108		1 SPL	Cup2	1						1	0,
109		1 SPL	Cup2	1						1	0,
110		1 SPL	Cup2	1						1	0,
901	baseline	1 RB	Cup1	1						1	0,
910	5 5 check	1 ICV	Cup1	1						1	0,
	0 blank	1 SPL	Cup2	1						1	0,
901	baseline	1 RB	Cup1	1						1	0,

Column 1: Cup Numbers

Positions 901-910 are for large tubes at the back of the sampler and are used for standards, checks and baselines (washes). The instrument can sample from these cups up to 35 times each.

There are 2 racks for the 2 ml sample cups. The one on the left (closest to the sampler wash cup) holds cup numbers 101-190. The second rack holds cup numbers 201-290. These cups should only be sampled one time each.

The instrument may have some sensitivity drift. Runs should not be longer than 90 samples. Larger sets of samples should be split into multiple runs.

Samples are loaded into the instrument in the pattern of 10 samples followed by a reagent blank, check, wash and a reagent blank. Repeat this pattern until the set is complete.

Column 2: Sample Names

Sample names should match the labels on the bottles to allow tracking of NEON samples from receipt to final data delivery.

The standards used in the calibration curve (anything labeled “STD(X)” in column 4) must be named exactly the same as they are named in the calibrant table in the method file.

To view the calibrant table, click the “edit method” button on the main toolbar, click file, open, and select appropriate method.

The range of the calibrant concentrations is chosen to cover the range of the samples. At least 5 calibrants are used with 7 or more being preferred.

Washes are DI water. Those with an “RB” in column 4 are used by the instrument software to account for baseline drift of the instrument. These should be run immediately before and after the calibration curve, at least every 10 samples during the run, and as the last sample of the run. Reagent blanks are either water or KCl depending on the analysis chosen.

Other washes are frequently run to clear any potential carryover after standards or samples with high concentrations. These will have “SPL” in column 4s or the software will not use them for baseline corrections.

Column 4: Type

SYNC-This is always the first cup in a run, it is a high standard (usually top standard) which will yield a large peak to let the software know the sample peaks have started.

STD- indicates a calibration standard to be included in the calibration curve. Do not use this for internal standards or checks unless you want them included in the calibration curve.

RB- indicates a baseline (wash). The instrument will use these peaks to correct for baseline drift so they should always be a 0 ppm standard.

SPL- unknown. Should be used for all samples and blanks.

ICV-indicates check standard

Other columns in the sample table

Column 3 - Number of replicate samples to be taken from the cup. Should be set to 1.

Column 5 - Cup size (Cup1 for 50ml tubes for the standards, Cup2 for samples in 2ml cups).

Column 6 - Dilution factor of the sample, usually 1 unless the sample has been diluted so that it will fall within the range of the standard curve. If the sample has been diluted, the dilution factor should be entered in this column (e.g. 1 part sample to 9 parts water would be a dilution factor of 10).

Columns 7-14 to remain unaltered to allow the CSV file to import into the software.

Quality Assurance

The range of the calibrant concentrations is chosen to cover the range of the samples. At least 5 calibrants are used with 7 or more being preferred.

Washes are DI water or KCl. Those with an "RB" in column 4 are used by the instrument software to account for baseline drift of the instrument. These should be run immediately before and after the calibration curve, at least every 10 samples during the run, and as the last sample of the run. For a nitrate run the reagent blank will be water whereas for the ammonium run the reagent blank will be KCl. The reagent blanks are taken from cup 901 which is also used as the zero for the standard curve.

Other washes are frequently run to clear any potential carryover after standards or samples with high concentrations. These will have a "SPL" in the type column so the software will not use them for baseline corrections.

Checks are standards prepared from a different source than the calibration standards. They are run after the calibration curve and about every 10 samples throughout the run. They are used to verify the calibration and to check for sensitivity drift of the instrument. Their concentrations should be between 20% and 100% of the highest concentration standard used for the run. Their sample name should include the nominal concentration.

NO₂ checks are standards prepared with sodium or potassium nitrite. They are run after the calibration to determine the efficiency of the cadmium reduction column. They should have nitrogen concentrations at least 50% of the highest calibration standard with 100% being the preferred value. Their label in the sample name column should include the nominal concentration.

Standard Operating Procedure: Processing NEON KCI Data

10/08/2018

Prepared by G. Beresford

EcoCore Analytical Services, Colorado State University

Exporting results from the Alpkem software

1. Open FlowView software and double click the instrument to get to the home screen.
2. Select View Results.
3. Select Open Results folder and choose file to be examined.
4. From the result window select export results and then the ... button beside the filename box.
5. Type desired filename and hit save.
6. Select the export button to move the data set to the Exported Files folder.
7. Close the results file and use a USB drive to copy the exported results to a networked computer.
8. Import the results into Excel. Use separate pages for the NH₄ and NO₃ results and label the pages "NH₄ raw" and "NO₃ raw". Leave these pages unaltered to serve as record of the results as they came from the instrument.

Processing the ammonium results

1. Copy the "NH₄ raw" page onto a new page labeled "NH₄ calc".
2. Move the lines containing the NO₃ calibrants, NO₃ checks, the NO₂ check, and the NO₃ baseline samples down to the bottom of the file.
3. Move the lines containing the NH₄ calibrants and the regression information down to about 30 lines below the last sample entry. Verify that the regression correlation is at least 0.995. *If it is not, the results are not valid and the samples will have to be rerun.*
4. Collect all of the lines containing the NH₄ check standards into a group below the sample results. Calculate a recovery for each check standard by dividing the measured value by the nominal value. The mean recovery should be between .95 and 1.05. *If not, the results are not valid and the samples will have to be rerun.*
5. Collect all of the baseline samples (type=RB) into a group below the check standards.
6. Check that all sample values are less than the highest standard. Flag any samples that are higher to be diluted and rerun.

Processing the nitrate results.

1. Copy the "NO₃ raw" page onto a new page labeled "NO₃ calc".
2. Move the lines containing the NH₄calibrants, NH₄ checks, and the NH₄ baseline samples down to the bottom of the file.

3. Move the lines containing the NO₃ calibrants and the regression information down to about 30 lines below the last sample entry. Verify that the regression correlation is at least 0.995. *If it is not, the results are not valid and the samples will have to be rerun.*
4. Move the NO₂ check result line below the sample results. Copy the line with the corresponding NO₃ calibrant result onto the line below the NO₂ check. Calculate the recovery by dividing the NO₃ calibrant result by the NO₂ check result. The recovery must be 0.9 or greater; *if not, the results are not valid and the samples will have to be rerun.*
5. Collect all of the lines containing the NO₃ check standards into a group below the sample results. Calculate a recovery for each check standard by dividing the measured value by the nominal value. The mean recovery should be between .95 and 1.05. *If not, the results are not valid and the samples will have to be rerun.*
6. Collect all of the baseline samples (type=RB) into a group below the check standards. Average the results. The average must be less than 0.02 or the MDL (whichever is greater). *If not, the results are not valid and the samples will have to be rerun.*
7. Check that all sample values are less than the highest standard. Flag any samples that are higher to be diluted and rerun.

Preparing data for submission to NEON

1. Copy sample names, measured results, and relevant run metadata for both ammonium and nitrate onto a summary page that matches exactly the NEON ingest datasheet template. *Do not include samples that have been flagged for exceeding the top standard as those must be re-run.*
2. Ensure all list-of-value fields are filled in using appropriate choices and any necessary quality flags are indicated using the QF fields.
3. Prepare an additional sheet with the batch QA (check standard) results for each run using the NEON batch QA ingest datasheet template. Ensure that run IDs match those in the sample data.
4. Save both the summary page and the batch QA page as separate .csv files and submit to NEON.