

**Laboratory of Medical Zoology Quality
Assurance/Quality Control Guidance for Arthropod-
borne Pathogens Testing***

Version 3.0 – April 8 2018

(*this document has been modeled after EPA QAQC guidelines
<https://www.epa.gov/sites/production/files/2015-07/documents/epa-qaqc-pcr.pdf>)

Table of Contents

1 INTRODUCTION

- 1.1 Code of ethics
- 1.2 Purpose
- 1.3 Scope

2 LABORATORY QUALITY ASSURANCE

- 2.1 Personnel
 - 2.1.1 Background and Training
 - 2.1.2 Outerwear
- 2.2 Laboratory background, design and workflow
 - 2.2.1 Laboratory background
 - 2.2.2 Laboratory design
 - 2.2.2.1 Reception Office (RO – 101 Fernald)
 - 2.2.2.2 In-Processing Terminal (I-PT – B1 Fernald)
 - 2.2.2.3 Nucleic Acids Extraction Lab (EXL – B1 Fernald)
 - 2.2.2.4 Reagent Prep Lab (RPL – B2 Fernald)
 - 2.2.2.5 Amplification and analysis room (PCRL – B2A Fernald)
 - 2.2.3 Workflow
- 2.3 Sample acceptance procedure
- 2.4 Equipment
 - 2.4.1 Real-time PCR Instruments
 - 2.4.2 Centrifuges
 - 2.4.3 Gel Electrophoresis Chambers
 - 2.4.4 Power Supplies
 - 2.4.5 Biological Safety Cabinets
 - 2.4.6 Pipettes
- 2.5 Disposables
 - 2.5.1 Pipette tips
 - 2.5.2 Sample and PCR tubes, Plates
 - 2.5.3 Gloves
- 2.6 Laboratory cleaning

3 REAGENTS, KITS, PRIMER SETS, AND ENZYMES

- 3.1 Reagents
- 3.2 Commercially available kits
- 3.3 Primer sets and hybridization probes
 - 3.3.1 Storage
- 3.4 Enzymes

4 ANALYTICAL PROCEDURES AND ASSESSMENT

- 4.1 Biological sample collection and processing
 - 4.1.1 Sample custody
- 4.2 Nucleic acid extraction
- 4.3 Polymerase chain reaction amplification
 - 4.3.1 Reaction volumes
 - 4.3.2 Primer and template concentrations
 - 4.3.3 PCR reagents and master mix preparation
 - 4.3.4 Selection of procedure parameters

- 4.4 Amplicon detection and confirmation
- 4.5 Method sensitivity, precision, and recovery
 - 4.5.1 Detection limits
 - 4.5.1.1 Detection limit of qPCR
 - 4.5.1.2 Detection limit of method
 - 4.5.1.3 Precision
 - 4.5.1.4 Recovery
- 5 INTERNAL QUALITY CONTROL PROCEDURES
 - 5.1 Positive controls
 - 5.1.1 Built-in Tick DNA Internal control
 - 5.1.2 PCR positive control
 - 5.1.3 Method positive control
 - 5.2 Negative controls
 - 5.2.1 Extraction negative control
 - 5.2.2 PCR negative control
 - 5.2.3 Method blank
 - 5.3 Reverse-transcriptase-qPCR quality control
 - 5.3.1 Contamination of RNA samples
 - 5.3.1.1 Sources of RNase contamination
 - 5.3.2 Preventing RNase contamination
 - 5.3.3 Quality-control samples for RT-qPCR analysis
 - 5.4 Corrective actions
 - 5.5 False-positive/false-negative prevention
 - 5.5.1 False positive prevention and detection
 - 5.5.1.1 Product carryover reduction
 - 5.5.2 False negative prevention
- 6 DATA RECORDING, RECORD KEEPING, AND DATA EVALUATION
 - 6.1 Data recording and record keeping
 - 6.1.1 Equipment
 - 6.1.2 Reagents, kits, primer sets, and enzymes
 - 6.1.3 Sample processing and analysis
 - 6.1.4 Quantitative real-time PCR

1 INTRODUCTION

Due to the ability to amplify small amounts of nucleic acid, PCR can be used to identify arthropod species and detect arthropod-borne pathogens that are difficult to culture *in vitro* or that cannot be cultured. However, the advantages of these techniques can be offset by the demanding assay protocols and the need to follow quality assurance/quality control (QA/QC) procedures carefully. These QA/QC procedures are necessary because the ability of PCR to produce many copies of target DNA creates the possibility of contamination by previously amplified products, which can lead to false-positive results. In addition, arthropod samples may contain PCR inhibitors, such as blood, which can lead to false-negative results. PCR protocols are standardized for analyses of arthropod samples, so it is essential to standardized QA/QC procedures.

The Laboratory of Medical Zoology at UMass-Amherst is specialized in the identification of arthropods, the detection and quantitation of arthropod-borne pathogens. The Laboratory of Medical Zoology provides a platform for multifaceted, multidisciplinary research at UMass-Amherst for the academic and public service.

The main goal of the Laboratory of Medical Zoology is to produce data that is scientifically valid, legally defensible, and of known and documented quality in accordance with nationally recognized standards.

1.1 Code of ethics

All members are to adhere to the following code of ethical practices:

All members will perform all work in a manner that merits full confidence and trust.

The Laboratory of Medical Zoology thereof will not engage in illegal practices, or cooperate with anyone so engaged.

The Laboratory of Medical Zoology thereof will ensure the integrity of their data by complete adherence to the laboratory QA/QC guidance, and will be diligent to expose and correct any errors that may be brought to light.

The Laboratory of Medical Zoology thereof will work and act in a strict spirit of honesty and fairness to clients, and in a spirit of personal helpfulness and fraternity toward fellow members.

The Laboratory of Medical Zoology thereof will not accept any of the following:

- Fabrication of data
- Misrepresentation of QC samples
- Non acceptable instrument calibration procedures
- Modification of samples to alter their characteristics
- Improper and unethical manual integrations
- Manipulation of analytical results
- Substitution of samples, files, or data

- Falsification of records or instrument readings
- Any other form of fraud

The Laboratory of Medical Zoology thereof will demonstrate the positive qualities of enthusiasm, diligence, responsibility, initiative, integrity, honesty, kindness, and patience in dealing with both clients and fellow laboratory members.

The Laboratory of Medical Zoology thereof, will advise clients of the probability of success before undertaking a project, and will not accept work that would constitute conflict of interest.

The Laboratory of Medical Zoology thereof will ensure the confidentiality of all data and information provided by the clients.

The Laboratory of Medical Zoology thereof will only perform testing services for which they are consistently demonstrated full compliance with high quality.

1.2 Purpose

The manual has been created as a resource for the Laboratory of Medical Zoology, QA/QC practices for arthropod sample collection, stabilization, storage, DNA and RNA extraction, qPCR data analysis, and for developing of new methods. Also it is a basis for principal investigators, researchers, veterinary practitioners, and students to evaluate the quality of PCR data for research projects and technical papers.

1.3 Scope

This manual is a general guidance for the Laboratory of Medical Zoology and method specific QA/QC procedures for real-time PCR data analysis of arthropod samples. This document does not address federal, state, and local regulations governing waste management, hazardous materials, and radioactive material. Furthermore, this guidance does not address related safety issues.

2 LABORATORY QUALITY ASSURANCE

This section provides guidance on general laboratory practices on day-to-day routine and provides recommendations for personnel, workflow, equipment, disposables, cleaning, the laboratory organization, and line of authority.

2.1 Personnel

The Laboratory of Medical Zoology director has the responsibility for the overall management and supervision of the laboratory and its personnel. He will interface with clients on all aspects of their projects including progress, problems, and recommended solutions. He will also work with the laboratory management and personnel in reviewing progress reports, analytical reports, financial reports, and QC reports. He is responsible for the review of all data generated in the laboratory for accuracy and interpretation. He is primarily responsible to ensure that all members have received the necessary level of training to make them capable of properly executing their duties.

The laboratory members are responsible for sample analysis, data processing and recording in accordance with the QA/QC guidance and established protocols. They are

responsible for calibration and preventive maintenance of instrumentation, data reduction, data review, and reporting of all out-of-control situations, as well as for initial corrective actions whenever necessary. Well documented training records are kept on file for each researcher in order to provide proof of proper training for each method they perform.

Personnel working in the laboratory performing qPCR analysis should meet background and training specifications outlined below and should follow the guidance provided concerning protective outerwear.

2.1.1 Background and Training

Technicians involved in arthropod testing should have undergraduate course work in biology, molecular biology, biotechnology, biochemistry, molecular genetics, or other course work that covers PCR and recombinant DNA theory and practice. However, commensurate job-related training and experience may be substituted. Hands-on training, including the review of standard operating procedures (SOPs) or manuals, should be completed for each technique under the supervision of experienced personnel (e.g., principal investigator, senior analyst, or manufacturers' representative).

New lab members receive orientation and skills training. New or established members may receive training on new methods given by the method developer.

Although the amount of time required for training will vary depending on the procedure and the technique, each member should demonstrate that they could successfully perform the assay through analyses of positive and negative control samples before being allowed to analyze arthropod samples without supervision. An initial demonstration of capability should be tested. Each member also is expected to be knowledgeable in laboratory safety and QA/QC procedures.

The Laboratory of Medical Zoology should maintain a training record for each researcher that documents the following:

- Dates and scope of TNA extraction/real-time PCR method training
- Initial demonstration of capability
- Proficiency test results for each analysis type
- Dates and scope of laboratory QA/QC training
- Dates and scope of laboratory safety training

2.1.2 Outerwear

UMass-Amherst provides laboratory coats program by Cintas. Coats are picked up weekly for laundering and returned to the lab the following week. Each lab member is designated two coats so that one is available in the lab while the second is out for laundering. If the primary coat becomes contaminated while the other is out for laundering, the primary coat can be decontaminated by either: a). autoclaving the current coat or, b). soaking it in a freshly prepared 10% bleach solution (bleach may damage labels), or c). have disposable lab coats available for "emergencies".

In addition to dedicated laboratory coats, powder-free gloves should be available in each laboratory room. Gloves should be discarded before leaving each room. Changing gloves reduces the possibility of contamination with template (the nucleic acid upon which the PCR is performed) or amplified nucleic acid. Gloves should be changed after working with seeded or environmental samples, after handling template or amplified nucleic acids, and after contact of the outside of the gloves with skin. The latter prevents introduction of enzymes prevalent on the skin, such as DNases and RNases, that degrade nucleic acids.

Guidelines for outwear and other necessary PPE are provided by UMass Environmental Health and Safety (<http://www.ehs.umass.edu/laboratory-and-research-safety>).

2.2 Laboratory background, design and workflow

2.2.1 Laboratory background

Laboratory facilities in the LMZ are open during normal business hours are Monday thru Friday, 9AM to 5PM. The LMZ Reception Office (101 Fernald) is open Monday thru Friday, 10AM-4PM on business days. The LMZ is closed on University holidays and/or cancellations (due to snow, for example).

The Laboratory of Medical Zoology (LMZ) laboratory and associated offices are located on three floors of Fernald Hall at the University of Massachusetts-Amherst.

LMZ Reception Office (RO; 250 ft²) is located on the first floor of Fernald Hall with seating for lab visitors and outreach consultation. All US Post (including arthropod packages) is received each morning as well as courier deliveries (FedEx and UPS) throughout the day in RO.

The In-Processing Terminal (I-PT; 2200 ft²) has designated bench space located in the basement laboratory (B1). This space has three stereoscope work stations equipped with cameras for photo documentation and imaging of specimens. Received specimens are matched with either online orders (TickReport™) or accompanying documentation (Agency Service), transferred to microcentrifuge tubes and entered into our internal database with unique identifiers cross referenced to contractor labels in the accompanying spreadsheets. LMZ codes are used throughout our workflow and final results are cross referenced with contractor info in final reports.

The large lab in the basement level comprises approximately 2200 sq ft. of laboratory space divided between three main rooms, B1 and B2a, B2b, with some smaller satellite rooms designed for particular procedures (cell/tissue culture; a darkroom; and microscopy room). Room B1 contains equipment and supplies utilized for processing of field and clinical specimens (biological safety cabinets, centrifuges, incubators, tissue lyser) in preparation for molecular biology experimentation. Archival storage (-20C, N=9 and -80C, N=1) are also partially located in B1.

Room B2 a is a large, well-equipped molecular biology laboratory designed to accommodate bench workstations for 8 bench scientists as well as common areas for gel and capillary, electrophoresis, conventional (N=4) and qPCR (N=4) instruments, bench top and refrigerated centrifuges (N=1). Two office spaces (~120 ft²/each) for the

lab manager and staff are adjacent to the B1 and B2 laboratories on the ground level of Fernald Hall.

B2 and B2a are equipped with a robotic liquid handling system (Eppendorf EpMotion) that allows for larger scale automated processing of samples as they are put through the process from clinical specimens to DNA extraction, PCR amplification, and downstream genotyping (SNPs, MSATs, nucleotide sequencing, etc.). All diagnostic assays are conducted in either 96- or 384-well format, providing high throughput while minimizing operator error. The lab also has three (3) Quantitative-PCR instruments (Agilent), a variety of electrophoresis units (vertical PAGE and DGGE apparatuses, horizontal agarose, CHEF, TGGE/DGGE tanks), water bath incubators, CO₂ incubators, benchtop and refrigerated centrifuges (N=5), refrigerators (N=4) and freezers (-20°C; N=9 and -80°C; N=1). Imaging equipment includes cameras for use with the labs stereo microscopes and the compound scopes (BF, DF, PC) and a gel documentation system. A gel-less capillary instrument eGENE HAD-GT12 capillary electrophoresis apparatus stains and detects bands in 96-well format in under 5 minutes to a resolution of 2 bp. Shared facilities within the department include access to scintillation counters, ultracentrifuges, and other common equipment.

PI, graduate student and staff offices are located on second floor of Fernald Hall in suite 209. Office space for students (~500 ft²), technicians and post-docs in an open common area with desks, carrels, conference table and book shelves. The PI's office (~300 ft² with kitchenette and conference area) is located in the corner of the suite. Space (~150 ft²) for senior research associates are located in this suite adjacent to the large common area. Networked printers are located in each of the main office areas. Office supplies, copier and printing materials are located in two satellite locations. At the center of the suite is a kitchenette/breakroom area, white boards/52" LCD screens and large (15') conference table with seating for 12.

A separate three story building (The Apiary) comprising ~2200 ft² is located nearby Fernald Hall and is currently in the LMZ inventory. This building contains offices on the top floor, a large open lab bench area on second floor with chemical hoods and a large walk-in growth chamber used for mosquito rearing, and a ground level area formerly used for vertebrate animals. LMZ has large 15 passenger van with all but one bench seat removed for conducting field work. Other vehicles are available through UMass fleet operations. All vehicles are maintained by physical plant personnel and gasoline is purchased at an on-campus supplier.

2.2.2 Laboratory design

The high sensitivity of PCR techniques requires that demanding assay conditions be followed. The Laboratory of Medical Zoology is designed and operated in a way that prevents contamination of reactions with amplified products from previous assays and cross-contamination between samples, both of which can lead to false-positive results.

Contamination between samples and from previous qPCR amplicons generated in the Laboratory of Medical Zoology could be significant potential source of invalid PCR results. Thus, the Laboratory of Medical Zoology has physically separate rooms:

- Reception Office (101 Fernald)
- In-Processing Terminal (B1 Fernald)
- Nucleic acids extraction room (B1 Fernald)
- Reagent preparation and QPCR set up room (B2 Fernald)
- Amplification and analysis room (B2a Fernald)

The workflow flows from top to bottom of this list to minimize risk of contamination. Equipment should not be moved between the rooms used for any PCR sample processing and analysis steps.

2.2.2.1 Reception Office (RO – 101 Fernald)

The reception office serves as the hub for public outreach and client engagement. It is also where all mail, including arthropod packages, is received each morning US Post and via courier deliveries throughout the day. Incoming specimens are passed to the In-processing terminal in the adjoining room.

2.2.2.2 In-Processing Terminal (I-PT – B1 Fernald)

Packages containing specimens are opened and inventoried in the in-processing terminal (I-PT). The IP-T has stereoscope work stations equipped with cameras for photo documentation and imaging of specimens. This room also contains an electronic label maker for labeling microcentrifuge tubes containing specimens that are carried on to the nucleic acid extraction room.

2.2.2.3 Nucleic Acids Extraction Lab (EXL – B1 Fernald)

The Nucleic acids extraction lab (EXL) is a space dedicated for organic extraction and purification of nucleic acids. The sample preparation room has dedicated adjustable pipettes and aerosol-barrier tips. Fresh gloves and laboratory coats are worn at all times. Gloves should be changed frequently including whenever it is necessary to leave the EXL area to control contamination from this room to other locations in the work flow.

2.2.2.4 Reagent Prep Lab (RPL – B2 Fernald)

The reagent preparation lab room contains a dedicated positive pressure hood for the preparation of PCR reagents, including master mixes and qPCR reactions set up. Addition of master mixes to PCR tubes should be performed in this hood in RPL. To prevent cross-contamination and to avoid repeated freezing and thawing, reagent-stock solutions should be aliquoted into 1ml “working” volumes and stored in -80C freezer. The extracted samples and extraction negative controls should be added to tubes/plates containing PCR master mix in this room. This room should have dedicated adjustable pipettes with plugged, aerosol-barrier tips. Clean gloves and laboratory coats should be worn at all times to control contamination from this room to any other location.


2.2.2.5 Amplification and analysis room (PCRL – B2A Fernald)

Two qPCR instruments are located in this room, which is designated for activities associated with PCR amplification and post-PCR analyses. Tubes or plates containing PCR reactions should be sealed (closed tubes or film sealed plates) while in this room. Gloves and laboratory coats should be worn at all times and removed before leaving the room to control amplicon contamination of other locations.

2.2.3 Workflow

We employ a unidirectional workflow to reduce the potential for contamination. Where appropriate color-coding is used to record progress of samples and instruments in the unidirectional workflow by designating colors specific to each laboratory room. In general, laboratory manuals and notebooks also should not be moved from room to room. Laboratories may want to consider using electronic data reporting to avoid this potential source of contamination.

Figure 1. Laboratory Workflow (Edits made to In-Processing + NA Extraction Sections of Figure 1)



2.3 Sample acceptance procedure

The Laboratory of Medical Zoology should assess arthropod samples when they are received at the Reception Office laboratory to verify that the sample was handled and preserved appropriately, the holding time requirement was met, and that all required sample collection information was recorded by the sample collector. After the sample is assessed, information on the date of sample receipt and sample condition should be recorded. The sample should be marked, logged, and tracked with a unique identifier.

2.4 Equipment

LMZ uses equipment suitable for the methods employed that function properly in order to generate reliable data. New instruments and equipment, or instruments and equipment that have undergone repair or maintenance, are to be calibrated before being used in casework analysis. To verify that equipment is functioning properly, the laboratory should have a schedule for maintaining equipment, which includes the setup, calibration, repair, record keeping, and normal operation of all equipment used in sample analysis.

2.4.1 Real-time PCR Instruments

Thermocyclers are essential to all PCR methods, and great care should be taken to

ensure that they are well-maintained and reliable. Real-time PCR instruments are equipped to perform fluorescence excitation and detection to monitor amplification throughout the PCR cycles. For this reason, the design is often different from the standard thermocycler, and calibration may be specific to the instrument design and should be carried out according to manufacturers' specifications. For example, Agilent MX3000P qPCR instruments eliminate the need for well-to-well signal correction by calibration or reference dyes. Uniformity of well-to-well measurements from the single-scanning fiber optic head design reduces the need for additional signal correction and calibration steps.

Qualification Test should be performed by using Agilent beta actin qualification plate every two years. LMZ should verify that the result cell of MX3000P reports a pass result when install a new qPCR machine and every two years.

2.4.2 Centrifuges

Separate centrifuges, including microfuges, should be used for pre- and post-PCR procedures. Centrifuges should have speed capabilities that meet the criteria required at different stages in the overall method. The manufacturers' instructions for calibration should be followed. The centrifuge should be balanced before use to increase bearing life and minimize vibrations that can unsettle concentrates.

2.4.3 Gel Electrophoresis Chambers

Chambers should be inspected before each use to ensure that electrodes and buffer tanks are intact, and that power supply electrodes fit snugly. Gel electrophoresis chambers should be rinsed several times with water after each use in the designated product room.

2.4.4 Power Supplies

Power supplies should provide voltage and current readings. Before use, electrodes should be checked to ensure that they fit snugly. Electrical cables, electrode wires, and cable connections should be checked for breaks, fraying, corrosion, or looseness, and replaced if needed.

2.4.5 Biological Safety Cabinets

PCR master mix preparation is conducted in a Class II biosafety cabinet equipped with UV lights for sterilization. This hood provides efficient laminar air flow to keep contaminants from entering the work area or the laboratory. Before use, hoods should be decontaminated using UV light for at least 15 minutes and cleaned with bleach or other effective nucleic acid inactivating agent. After a spill or accident, hoods also should be cleaned with the nucleic acid inactivating agent and then decontaminated further using UV light. Airflow and HEPA filtration in all hoods should be monitored and certified per manufacturers' recommendations annually.

2.4.6 Pipettes

Due to the numerous small volume transfers involved in PCR methods, automatic pipettes, fixed-volume, adjustable, positive-displacement, and/or micropipettes should

be used, and should be calibrated annually by the manufacturer or a technician using the tips commonly used in the laboratory. Documentation on the precision of the recalibrated micropipettes is retained and kept on file in the laboratory's QC logbook. Each pipette should be sterilized according to manufacturers' recommendation on a regular basis or whenever contamination is suspected.

Pipettes failing to meet calibration recommendations will be decommissioned and/or sent to manufacturer for repair or replacement.

2.5 Disposables

Disposable materials used in PCR analysis include pipette tips, sample tubes, PCR tubes, and gloves. To reduce the contamination and degradation of the target nucleic acids, disposable materials should meet the standards discussed in sections 2.5.1 to 2.5.3.

2.5.1 Pipette tips

Special tips for PCR analysis include barrier tips and aerosol-resistant tips, both of which minimize cross-contamination of samples during pipetting. Pipette tips for PCR analyses should be lot-certified, RNase-free, DNase-free, and pyrogen-free.

2.5.2 Sample and PCR tubes, Plates

LMZ uses polypropylene tubes and plates that are lot-certified DNase-, RNase-, and pyrogen-free. Thin-walled tubes provide the best heat transfer, ensuring that the reaction volume reaches its specified temperature in the shortest amount of time, thereby improving specificity and reproducibility. Tubes and plates containing stored samples and reagents should be centrifuged briefly before opening to ensure that all liquids are at the bottom.

2.5.3 Gloves

Disposable gloves should be available in each section of laboratory. Gloves should be changed before leaving and entering each section of the laboratory and each time contaminating DNA is potentially encountered. In addition to reducing potential contamination from samples, wearing gloves may protect the researcher from potential chemical exposure and prevent sample contamination due to human DNase and/or RNase.

2.6 Laboratory cleaning

All work surfaces should be cleaned after each use with 0.6% sodium hypochlorite (NaOCl). The NaOCl solution should be prepared fresh daily by diluting commercial bleach 1:10 in water and adjusting the pH to 7. This solution inactivates pathogenic agents and destroys nucleic acids. Residual bleach, which may cause pitting in stainless steel counter tops and hoods, may be removed with a 70% ethanol rinse or equivalent. Commercial products that are specifically designed for removing nucleic acids and nucleases also can be used for surface cleaning. Thermocyclers and centrifuges should be cleaned with the diluted bleach solution whenever contamination is suspected. Pipettes should be cleaned according to manufacturers' instructions.

3 REAGENTS, KITS, PRIMER SETS, AND ENZYMES

3.1 Reagents

The reagents used in PCR amplification can be purchased, and should be molecular grade. These reagents should be stored according to the manufacturers' recommendations. All reagents from new lots should be tested to ensure that they work properly by running a PCR positive control. Care should be taken to ensure the reagents are maintained contamination-free. All reagents should be clearly labeled with name, expiration date, and relevant safety information. Reagents from different lot numbers should not be interchanged without prior functional validation.

Molecular-grade water or its equivalent from commercial sources (DNase/RNase-free) should be used for all assays. Diethylpyrocarbonate (DEPC) treatment can be used to eliminate RNase from water used in RNA analysis. Reagent water is treated with a solution of 0.1% DEPC for several hours and then autoclaved to degrade the DEPC completely. Proper autoclaving is necessary, because trace amounts of DEPC in a solution will lead to the modification of the purine residues in RNA by carboxymethylation. This leads to downstream effects in RNA experimentation.

3.2 Commercially available kits

Many types of commercial kits are available for isolation of DNA and RNA and the purification of nucleic acids to remove contaminants. A copy of the manufacturers' specifications and procedures should be kept in the laboratory's binder. The effectiveness of each new kit should be evaluated before use for the analysis of field samples by running an appropriate positive control.

3.3 Primer sets and hybridization probes

Primer sets should be designed specifically to prime the amplification of a portion of a target nucleic acid of interest. Hybridization probes are specific for nucleotide sequences that are internal to the sequences of the primers and which are used to confirm the amplification of the target. As a general rule, primers are characterized by the following:

- Length of 18 to 27 base pairs
- No homology within or between primers, especially at the 3' end to avoid primer-dimer formation
- No guanine-cytosine (GC) stretches greater than two base pairs at the last five nucleotides at the 3' end
- GC content of 30% to 80%
- Melting temperatures (T_m) of 58-60 °C, and probes with a T_m value of 10 °C higher
- No hairpin loops with an energy of -0.5 kcal/mol or less
- It is preferable to have primers spanning exon-exon junctions in the cDNA sequence

- Maximum amplicon size should not exceed 400 bp (ideally 50-150 bases)

The probes should not have runs of identical nucleotides (especially four or more consecutive Gs), there should be more Cs than Gs, and not a G at the 5' end. The choice of probe should be made first.

The Laboratory of Medical Zoology uses software for primer design. New primers and probes should always be tested experimentally. The specificity of a chosen sequence should be evaluated using BLAST. Specificity testing should be conducted as described for the PCR positive control by substituting the heterologous sample in place of the target organism or nucleic acid. Primers and/or probes should be redesigned if any heterologous sample is positive. Primers and probes should be selected based on the best theoretical considerations, and negative features should be documented, and results interpreted accordingly.

Every lot of new oligos should be checked for contamination by being used in a PCR negative control. No positive results should be found. Primers and probes should be added to the PCR master mix in the reagent preparation area.

Functional validation also should be performed on new lots of primers and probes by comparing their performance against older sets of known quality.

The PCR efficiency (E) is calculated from the slope of a standard curve (S) using the formula, $E=10^{-1/\text{slope}}$. The samples used for generating standard curves should be carefully prepared dilutions of a stock solution containing a known number of purified PCR template nucleic acids or total nucleic acids extracted from the target organism. A

PCR efficiency of 1 (or 100%) is obtained when each target sequence present in the PCR reaction doubles during each round of amplification. A high efficiency is needed to obtain accurate and reproducible results. A low efficiency indicates that primers and/or probes are performing poorly. No single efficiency (as a percentage) has been set as a limit for data to be acceptable, although 80-120% is an accepted norm.

The default concentration of primers and probes is 200 uM in the Laboratory of Medical Zoology.

3.3.1 Storage

Most oligos and DNA templates should be stored at -20°C or -80°C in either TE buffer (10 mM Tris-HCl and 0.1mM EDTA, pH 8.0) or molecular grade water. PCR products may also be stored at -20°C or -80°C. RNA templates should be aliquoted and stored at -80°C. The pure, concentrated oligos should be stored in the original tube from the manufacturer and labeled with the primer name and concentration. To minimize the chance of contamination and degradation, these concentrated stocks should not be used on a regular basis. Diluted working stocks should be made for each oligo, and these working stocks should be used for all experiments. The sensitivity of primers can be checked by running a PCR positive control using the old and new primers. If degradation has occurred, new primers and probes should be prepared.

3.4 Enzymes

Enzymes should be purchased from a commercial source to ensure purity. The

Laboratory of Medical Zoology selects vendors that provide QA information with the enzymes. After receipt from commercial sources, each new lot of enzyme should be compared with old lots using known controls. The new enzymes should be rejected or the concentration adjusted, if performance is significantly below that of the proven lot. The researcher should have only one container open at a time when working with enzymes.

Pipette tips should be discarded after each use to prevent cross-contamination. Enzymes typically are stored at -20°C, and should never be left at room temperature in the laboratory. Insulated bench-top coolers or ice can be used to keep the enzyme cold in the laboratory, when used on the bench top.

4 ANALYTICAL PROCEDURES AND ASSESSMENT

The Laboratory of Medical Zoology maintains standard operating procedures that accurately reflect all laboratory activities such as assessing data integrity, corrective actions, handling customer complaints, and all test methods. When developing a new PCR system, the laboratory should design and select the individual components of the system, to optimize its ability to recover and detect the target genes. The laboratory should consider how samples will be collected and processed for the system, and assess the performance of the individual system components and the entire analytical process before using the system for biological sample analyses.

The laboratory shall have and follow written analytical procedures, and have a standard operating protocol for each analytical technique used. The procedures shall include reagents, sample preparation, extraction, equipment, and controls that are standard for DNA/RNA analysis and data interpretation.

4.1 Biological sample collection and processing

Sample collection and transport conditions impact the results of PCR analysis. So, The Laboratory of Medical Zoology should advise principal investigators, researchers, veterinary practitioners, and students to the collection, storage, and transportation of samples to be analyzed. Acceptable sample volumes, sample handling, and the time samples can be held before being processed should be described, in order to preserve the integrity of the target nucleic acid sequence. The sample processing is performed to extract the nucleic acid of interest, while minimizing co-purification of potential contaminants and components of the matrix that may inhibit PCR. The extraction procedure also should stabilize the target nucleic acid from nucleases and reduce the sample to a volume small enough to be analyzed by PCR. Positive controls that demonstrate that the target organism or nucleic acid has been successfully isolated and that the PCR inhibitors have been removed.

4.1.1 Sample custody

The Laboratory of Medical Zoology has a sample acceptance policy that outlines the circumstances under which samples will be accepted. The policy requires complete documentation, including the sample identification, the location, date and time of collection, collector's name, preservation type, sample type and source, and any special remarks concerning the sample.

If a sample discrepancy, such as a broken or missing sample is observed at check-in, a statement to that effect is written into the data base remarks section. The client must be notified and approve of any changes made. At this time, the laboratory manager is notified so that the problem can be addressed.

All samples received are recorded into the data base with the following information: client name, client address, analytical parameter requested, and laboratory work order number. The laboratory work order number is a sequential number that is unique to each sample. Samples are processed through the laboratory by their unique laboratory work order number. Special instructions about the samples are written and are entered into the data base in the comments section.

Samples which specifically require thermal preservation are stored under refrigeration. Samples are stored in a manner that prevents cross contamination.

All samples, digestates and extracts or other sample preparation products are disposed of in accordance with Federal and State laws and regulations.

4.2 Nucleic acid extraction

Some methods require the lysis of the biological samples and the isolation of its nucleic acid before proceeding to PCR. Extracted nucleic acids should be assessed as part of the entire analytical process. The efficiency of nucleic acid extraction varies with sample type and extraction procedure. Ideally, an isolation technique should meet the following goals:

- Demonstrate a high efficiency of target nucleic acid recovery
- Maintain nucleic acid integrity and minimize fragmentation
- Provide sufficiently pure nucleic acid, free from PCR inhibitors
- Minimize the use of dangerous chemicals
- Be repeatable

4.3 Polymerase chain reaction amplification

To analyze biological samples using PCR, the first step is to select the type of PCR that is most appropriate for the analysis being conducted. The performance of the PCR amplification step should be assessed separately first, then as part of the entire analytical method. Aspects that should be considered in the selection of the appropriate technique include PCR type, enzyme type, primer and probe design, and reaction parameters, such as thermal cycling temperatures and reagent concentrations.

The Laboratory of Medical Zoology employs following PCR:

- Conventional PCR uses a thermostable DNA polymerase to amplify a region of target DNA defined at each end by a specific primer. The exponential replication of the same target sequence produces enough amplicons for use in subsequent analyses.
- qPCR detects and measures the amplification of target nucleic acids as they are produced. qPCR requires the use of primers and an oligonucleotide probe

labeled with fluorescent dyes or an alternative fluorescent detection chemistry.

- A computer is used to monitor the fluorescence increase and to calculate a quantification cycle (C_q) value. qPCR allows quantification of template by absolute standard curve method, relative standard method, or comparative C_q method. There are several different fluorescent detection chemistries used for real-time PCR, including the following:
 - Taqman probes
 - SYBR Green
- Reverse transcriptase-qPCR is used to amplify RNA target sequences, such as messenger RNA and viral RNA genomes. This type of PCR involves an initial incubation of the biological sample or control RNA with a reverse transcriptase enzyme and a DNA primer. DNA primers that are used commonly include oligos dT (an oligos consisting of only thymidine residues), random hexamers (primers made of six random nucleotides), or a specific primer.
- Nested PCR or semi-nested PCR is used to increase both the specificity and sensitivity of target sequence detection. The use of a second amplification step with the "nested" primer set results in a reduced background and increased amount of amplicon. However, there is an increased chance of carryover or cross-contamination when taking product from the first round of PCR and putting it in the tubes for the second round of PCR. Additional steps and precautions may need to be taken to reduce the chance of sample contamination and false-positives, including the following:
 - Never opening more than one tube at a time
 - Using a separate thermocycler
 - Adding additional negative control

4.3.1 Reaction volumes

PCR volumes typically range from 10 to 100 μ L and may be dependent on the type of thermal cycling instrument and method being used. Increasing the volume of sample added to a reaction may increase the probability of a positive result at low target concentrations because more samples can be tested per reaction, but may also increase inhibition. The default reaction volumes is 16 μ L in the Laboratory of Medical Zoology.

4.3.2 Primer and template concentrations

The primer concentrations used in each newly developed PCR assay should be optimized to obtain maximum amplification efficiency. The addition of too much DNA (either template or nontemplate) to the reaction may be inhibitory.

4.3.3 PCR reagents and master mix preparation

PCR kits containing all reagents, except primers and templates, are purchased from commercial sources and used by the Laboratory of Medical Zoology after being

optimized for the particular assay. In addition to required reagents, enhancers that act on different aspects of the reaction may be added to the PCR master mix, and these include dimethyl sulfoxide (DMSO), betaine, tetramethylammonium chloride (TMAC), and bovine serum albumin (BSA). Other proprietary enhancers and additives are available from various manufactures. Each additive should be empirically tested with each system.

A master mix containing optimized reagents should be prepared in the reagent preparation room. During master mix preparation, pipette tips should be changed after handling each reagent. The master mix should be prepared on ice, unless the manufacturers' instructions specifically state that the reagents can be mixed at room temperature. The master mix should be aliquoted into reaction tubes in the reagent preparation room. The laboratory can prepare a large quantity of master mix and store aliquots at -80°C after demonstrating that the freeze/thawing of aliquots does not reduce sensitivity.

4.3.4 Selection of procedure parameters

The exact parameters of the PCR should be selected based on the target nucleic acids that are to be amplified. Parameters that need to be evaluated for successful PCR amplification include thermocycling conditions, reaction volumes, template concentration and the concentration of PCR reagents, including that of primers. After optimal conditions are established, they should be documented for the method and followed thereafter.

4.4 Amplicon detection and confirmation

PCR results must be confirmed. Probe-based qPCR is a method employed by the laboratory, which confirms the presence of a specific product by detection of the fluorescence dye.

Note: Due to the small size of some amplicons in qPCR detections, confirmation by sequencing may not be possible.

4.5 Method sensitivity, precision, and recovery

After the individual components of the method have been optimized, the entire method should be evaluated to determine the method's detection limits (sensitivity), precision, and recovery.

4.5.1 Detection limits

Target organisms are often present at low numbers, and some of these may be able to produce infection at low concentrations. Also, expression of selected genes while monitoring their changes can be at very low levels. Therefore, an acceptable PCR method should be applied to detect a level of organisms or genes targeted that is biologically significant. The level may be based upon acceptable risk levels, predicted environmental concentrations, or concentrations that produce an infectious dose. In order to know whether an acceptable level can be detected, a detection limit must be determined for each system.

The systems employed by the LMZ should have the detection limit of at least 5pg of

DNA and RNA gBlock targets with the confidence level typically from 95% to 99%. For presence/absence methods, the detection limit is the minimum concentration of a target that produces a positive response with a given level of confidence. There are many uncertainties that can affect the detection limit. Some of them are:

- The type of target nucleic acid being detected (e.g., DNA, mRNA, tRNA, rRNA, etc.)
- The secondary structure and the GC content of the nucleic acid target molecule
- The matrix from which the organism is isolated
- The detection of microbes that are inactivated by physical and chemical disinfectants

The detection limits can be determined by analysis of replicate, seeded samples containing increasingly lower levels of the target organism or target nucleic acid. The number of replicates that should be analyzed at each dilution level will depend on the variability of the procedure, the accuracy and variability of the spike enumerations, the desired level of confidence in the detection limit, and the application of the method.

Errors associated with serial dilutions should be minimized by using appropriate volumes that can be accurately measured and by thoroughly homogenizing stocks immediately before removing an aliquot for the next dilution. Multiple replicates of each dilution should be enumerated independently of the PCR analysis to assess the accuracy of the dilutions.

4.5.1.1 Detection limit of qPCR

The PCR detection limit reflects the sensitivity of the PCR procedure, which includes the sensitivity of the primers and probes as well as the preparation of the master mix and the optimization of thermocycling conditions. It does not include reductions in sensitivity by any procedures used to concentrate organisms in the sample and isolate nucleic acids. For this reason should be determined on the basis of nucleic acid concentrations rather than organism concentrations.

4.5.1.2 Detection limit of method

The method detection limit is the lowest amount of target that can be reproducibly detected. It reflects not only the sensitivity of the qPCR, but the efficiency of the procedures used to detect target from the biological sample and the efficiency of the procedures for recovering nucleic acid. Thus, a method's detection limit will depend on a range of factors, including the efficiency of the collection, concentration, lysis of sample, and extraction of nucleic acid from the biological sample, in addition to the detection limit of the PCR itself.

The method detection limit can be determined by processing standardized matrix, or biological samples that have been spiked with different concentrations of the target organism prior to any sample processing.

4.5.1.3 Precision

Good PCR methods will have a high precision. Precision is a measure of how closely values from replicate measurements of a sample agree with each other. For quantitative PCR analyses on biological samples, method precision can be expressed as the relative standard deviation (RSD) or the relative percent difference (RPD). The RSD is the standard deviation from three or more replicate samples divided by the mean and multiplied by 100%. The RPD is the absolute difference of duplicate samples divided by the mean and multiplied by 100%. The RSD or RPD decreases with increasing precision and increases with increased variability. Variability (i.e., low precision) may result from differences in performance, reagents or equipment, or sample characteristics.

4.5.1.4 Recovery

In addition to high precision, good analytical methods will have low bias. Bias is a measure of disagreement between the concentration of a DNA or RNA template as measured by a method and the true concentration in the sample that was tested. A PCR method's recovery of the target sequence of interest is an important factor for characterizing method bias at levels higher than the detection limit. Bias can occur at any component of a PCR method. If a method is designed to detect a group of related organisms, members of the group may have different recoveries and thus different biases through the sampling protocol, nucleic acid extraction, and the PCR steps. Bias at the PCR level can occur from the favoring of the amplification of certain target sequences due to properties of the target, the flanking sequences, or the overall genome. PCR bias may be seen for mixed-template PCRs, including reactions using internal controls, and may make interpretation of results difficult, because the final results do not reflect the original makeup of the target templates. Bias usually results when the amplification efficiencies of the target sequences are not the same. Target sequences may amplify at different efficiencies because the sequences are not equally accessible to primer hybridization after denaturation, primer-template hybrids do not form with equal efficiency for all templates, and/or the polymerase acts on the templates with different efficiencies. One factor that may effect primer-template hybridization is the GC content of the template. PCR bias may be reduced by running a low number of PCR cycles and by using high template concentrations, when possible.

Recovery of the PCR analysis and the entire method can be assessed by the analysis of seeded samples. For qPCR analyses, recovery is determined as the total amount of the template found in a spiked sample minus the background (i.e., amount detected in the non spiked sample) divided by the amount of the template spiked into the sample. To determine recovery for presence/absence methods, serial dilutions of the spiked sample should be performed and five replicates at each dilution analyzed. The amount of template in the spiked and non spiked samples should then be calculated.

Results that are obtained using methods that have low recoveries should be qualified as potentially biased. Also the recovery of target organisms that are endogenous to a matrix (e.g., enteric pathogens associated with fecal matters) may not be the same as that measured by spiked samples (e.g., where the organisms are not associated with other material). Unless it can be shown that the recoveries of endogenous and spiked organisms are similar, sample results should be recorded as biased.

5 INTERNAL QUALITY CONTROL PROCEDURES

The Laboratory of Medical Zoology employs quality control samples to assess the validity of the analytical results of all samples. Determination of the validity of all sample results is based on the acceptance criteria being met by the control sample. The control samples are analyzed in the same manner as the samples. They are interspersed with the samples at frequencies that are specified by the appropriate standard operating procedure. The frequencies may be altered to comply with client-specific requirements. The data acquired from QC procedures is used to estimate the quality of analytical data, to determine the need for corrective action in response to identified deficiencies, and to interpret results after corrective action procedures are implemented.

The laboratory should analyze positive and negative QC samples on a routine basis to demonstrate adequate performance of PCR based methods.

5.1 Positive controls

Positive controls are analyzed to verify that the method is capable of adequately recovering and amplifying the target. The concentration of the sequence of interest in these positive controls should be 10 to 100 times higher than the defined detection limit of the PCR. A positive control is considered to be acceptable if the DNA/cDNA of interest was amplified by PCR. Precautions should be taken to avoid contamination of biological samples with the positive control template. Positive control preparation should be physically separated from biological sample and negative control preparation, and positive control samples should be handled last.

A summary of the recommended positive QC samples is as follow:

Description	Purpose	Frequency
Built-in Tick DNA Internal control	Verify that the sample DNA/RNA quality	Each sample
PCR positive control	Verify that the PCR master mix and reagents were prepared correctly to produce amplification of the target nucleic acid	One or more per PCR batch
Method positive control	Verify that the entire method is performing properly and is capable of amplifying the target nucleic acid from the organism of interest	With every sample batch

5.1.1 Built-in Tick DNA Internal control

We have a built-in tick DNA internal control for all hard ticks submitted, a built-in tick RNA control for *Ixodes scapularis* and *Ixodes pacificus*, a built-in tick RNA control for

Amblyomma americanum and a built-in tick RNA control for *Dermacentor variabilis*. Built-in Tick DNA Internal control is used to verify that the sample DNA/RNA quality is good enough for PCR reaction. It is also used to verify that interfering constituents from an environmental matrix carried over from the isolation of the organism or nucleic acids do not inhibit the PCR.

5.1.2 PCR positive control

PCR positive controls are used to verify that the PCR master mix and reagents were prepared correctly in order to produce amplification of the target nucleic acid. This type of positive control is run with each PCR batch. A PCR batch is defined as a group of samples that are processed and amplified at the same time under the same conditions, using the same PCR master mix, and in the same thermocycler. PCR positive controls are prepared by the addition of an exogenous control to the master mix. Exogenous controls can be any of the following:

- A purified total nucleic acid extract from the organism or the biological sample containing the sequence of interest.
- The whole organism, which may be used when the nucleic acid target of interest can be released from the seeded organism by heating before or during PCR.
- A specific nucleic acid fragment containing the entire sequence to be amplified, including primer binding sites.
- A cloned DNA fragment containing a modified form of the target sequence.
- A heterologous sequence that has been previously shown to be amplified with efficiency that is comparable to the target sequence.
- RNA transcribed from a cloned DNA fragment containing the target sequence or a modified target sequence (for use in RT-qPCR applications).

For quantitative PCR methods, the PCR positive control is evaluated by determining the total amount of the target nucleic acid or organism in the control divided by the amount added to the reaction as a spike.

5.1.3 Method positive control

The method positive control is used to verify that the entire method is performing properly. This control should be performed by analyzing a reagent biological sample spiked with known quantities of the target organisms prior the start of sample processing. For quantitative PCR methods, recovery of the method positive control is determined as the total amount of the analyte found in the sample divided by the amount of the control analyte added into the sample as a spike. The recovery should be expressed in terms of the units used to measure the concentration of the spike. At a minimum, method positive controls should be performed for each sample batch. A sample batch is defined as a set of test samples setup and processed together through all steps of the method leading to PCR.

5.2 Negative controls

Negative controls using each primer set and a probe should be analyzed to verify that no contaminating nucleic acid has been introduced into the master mix or into samples during sample processing. These negative controls are considered acceptable if no amplification of nucleic acids is detected. The negative controls are as follow:

Description	Purpose	Frequency
Extraction negative control	Verify that no contaminating nucleic acid has been introduced into the DNA/RNA extraction	1 negative per extraction set (typically 29 samples)
PCR negative control	Verify that no contaminating nucleic acid has been introduced into the master mix	At least 1-5% of the number of field samples analyzed per PCR batch
Method blank	Verify that no contamination has been introduced throughout the entire sample processing	With every sample batch

5.2.1 Extraction negative control

Nucleic acid extraction negative controls are used to verify that no contaminating nucleic acid has been introduced into DNA/RNA extraction. These controls are prepared during nucleic acid extraction step. They are prepared as separate samples to which aliquots of molecular-grade water or buffer are added to the Nucleic acid extraction step in place of arthropod sample. The number of nucleic acid extraction negative controls should be at least 1 per 29 arthropod samples.

5.2.2 PCR negative control

PCR negative controls are used to verify that no contaminating nucleic acid has been introduced into the master mix. These controls are prepared when template is added to the master mix. They are prepared as separate samples to which aliquots of molecular-grade water or buffer are added to the master mix in place of target nucleic acid or sample. A negative result with this control indicates that the master mix and final processing reagents are not contaminated. The number of PCR negative controls should be 1-5% of the field samples analyzed per primer set per PCR batch. A PCR batch is defined as a group of samples set-up and amplified at the same time under the same conditions, using the same PCR master mix, and in the same thermocycler.

5.2.3 Method blank

The method blank is designed to check for contamination throughout sample processing and PCR analysis. This control is performed on a sterile reagent water sample that is processed with the test samples using the same preparation, extraction, sample transfer, and PCR procedures as the test samples. At a minimum, method blank samples should be performed once per batch. A sample batch is defined as a set

of test samples processed together through all steps of the method leading to PCR. Method blank samples may need to be run more frequently, depending on the data quality needs or method and monitoring program requirements.

5.3 Reverse-transcriptase-qPCR quality control

Reverse transcriptase PCR reaction is a process in which single-stranded RNA is reverse transcribed into complementary DNA (cDNA) by using total cellular RNA or poly(A) RNA, a reverse transcriptase enzyme, a primer, dNTPs and an RNase inhibitor. The resulting cDNA can be used in RT-qPCR reaction. Three types of primers can be used for RT reaction: oligo (dT) primers, random (hexamer) primers and gene specific primers.

5.3.1 Contamination of RNA samples

Working with RNA, prevention, detection and elimination of nuclease contamination, particularly ribonuclease (RNase) contamination, should be constant. RNase plays an important role in nucleic acid metabolism, and is found in both prokaryotes and eukaryotes, and in practically every cell type. RNases tend to regain their native structure and partial function after being cooled to room temperature in the absence of a denaturant, and consequently can retain activity after freeze-thaw cycles and even after autoclaving. The robust nature of these enzymes makes them refractory to many methods of decontamination.

5.3.1.1 Sources of RNase contamination

RNA extracted from biological samples can be contaminated by RNase from several sources:

- Bodily fluids such as perspiration
- Tips & tubes (always use tips and tubes that have been tested and certified RNase-free)
- Water and buffers
- Laboratory surfaces
- Endogenous RNases
- RNA samples
- RNA storage
- Chemical nucleases
- Enzymes

5.3.2 Preventing RNase contamination

Some basic precautions need to be taken when working with RNA. These include:

- Wearing gloves throughout experiments and changing them after touching skin, door knobs, and common surfaces

- Having a dedicated set of pipettors that are used solely for RNA work
- Using tips and tubes that are tested and guaranteed to be RNase-free
- Using RNase-free chemicals and reagents
- Designating a "low-traffic" area as an "RNase-free zone"

Eliminating contamination from laboratory surfaces, glassware, metal ware, and pipettors by wiping it down with a commercially available RNase decontamination solution such as Ambion's RNaseZap, 10% solution of bleach, or NaOH. Spray or wipe the inside of biological safety cabinet or clean area used for sample processing with RNaseZap, or 10% bleach. In hoods, turn on the UV light and wait about 15 minutes before proceeding.

5.3.3 Quality-control samples for RT-qPCR analysis

It is widely known that extracted RNA can subject of degradation in various samples. RNA degradation especially taking place in paraffin-embedded tissues, but degraded RNA to various degrees has been reported in freshly collected biological samples. However, when using state of the art technologies, such as immediate transfer of biological samples into a solution that inhibits RNases and stabilizes RNA, degraded RNA is usually, but not always, absent after RNA isolation. To correct false gene expression measurements in potentially degraded samples, it is necessary to quantify the height of degradation. The Laboratory of Medical Zoology uses an adequate house-keeping gene with almost constant copy number in tissues of different origin.

The RT-qPCR is a complex assay and all physical and chemical components of the reaction are interdependent. They must be considered carefully when optimizing the specificity, sensitivity, reproducibility or fidelity of the reaction. To assess every step following controls should be employ:

SAMPLE TYPE	DESCRIPTION	PURPOSE
Negative-process control	Sterile water analyzed with each batch of samples beginning at inhibitor removal	Cross-contamination of samples during inhibitor removal, RT-qPCR
Spiked negative-process control (SPUD assay)	Total RNA added to a duplicate negative process control with each batch of samples before RT-qPCR	Verify proper completion of RT-qPCR reaction
Matrix spikes	Total RNA added to duplicate biological samples before RT-qPCR	Indicates whether inhibitors in environmental samples were present
PCR positive controls	RT-qPCR products from positive control added to second step reaction	Verify the proper completion of the RT-qPCR
PCR negative controls	Sterile water added at the RT-qPCR step	Contamination during RT-qPCR

5.4 Corrective actions

If any positive control failures occur, all samples associated with the control should be considered invalid, and negative field samples should be listed as potentially false-negative samples. If amplification of a positive control fails to produce the specific amplification product, the integrity of the control and the PCR design should be examined to determine the reason for the failure. When determined, the reason for the failure should be documented and the controls and samples re-run. Samples associated with matrix spike or inhibition positive controls that have low or no recovery should be qualified as potentially biased or false negatives, as appropriate, and the matrix should be evaluated to determine the cause of the interference, if possible.

For some samples, it may be possible to correct for bias in test results based upon the recovery as determined by the matrix spike. Corrections cannot be made if the PCR inhibition positive control shows evidence of inhibition. Corrections also cannot be extrapolated to samples for which matrix spike data are not available.

If PCR negative controls or method blanks produce specific amplification products, all samples associated with the failed controls should be considered invalid, and all positive samples should be listed as potentially false-positive samples. The source of contamination should be identified and eliminated. If the source of the contamination cannot be identified, additional types of negative controls should be added at various steps in the method to determine where the contamination is being introduced. Three such controls are as follow:

- Equipment blank. This control is performed by passing sterile reagent grade water or buffer through the equipment and processing the water as if it were a PCR negative control. If this control is found to be positive, all analysis should cease until the source of the problem is identified. Equipment blanks should then be run more frequently until it is shown that the problem has been corrected.
- Wipe test. This control is performed by wiping an area, pipette, or other equipment with a sterile, nucleic acid-free, gauze, then suspending the wipe in molecular-grade water and processing it as a PCR negative control. To determine whether nucleic acids detected are part of the normal background or the source of sample contamination, wipe tests should be performed routinely to establish a baseline for comparison.
- Room QC controls. This control is prepared by adding sterile water in place of template to a master mix in the reagent preparation room. These tubes are placed in working areas in the reagent preparation and sample preparation rooms and opened for 15 minutes. The tubes are then closed and run as a PCR negative control. These controls may be run on a monthly basis if contamination is a concern. If any of these controls are positive, all work should cease until the source of the contamination can be identified and corrected.

5.5 False-positive/false-negative prevention

Specific procedures should be implemented to prevent false-positive and false-

negative reactions from the techniques used. Details on preventing false positives and false negatives are provided in sections 5.5.1 and 5.5.2.

5.5.1 False positive prevention and detection

False positives can occur from contamination and non-specific amplification. Contamination introduced by personnel, positive controls, or positive samples (i.e., cross-sample contamination) can be minimized by the guidance provided in sections 2 and 3, and especially by physically separating work areas.

5.5.1.1 Product carryover reduction

Product carryover (contamination with DNA that has been previously isolated and/or amplified) can result in contamination of reagents or samples. Several approaches are available to prevent false-positive reactions through carryover contamination including the following:

- Enzymatic inactivation
- Photochemical inactivation (addition of psoralen after the PCR links the double-stranded DNA and prevents amplification in subsequent PCR)
- Hydroxylamine treatment (hydroxylamine is added after the PCR and prevents cytosine from bonding with guanine in subsequent PCR)

The most widely used approach, enzymatic inactivation, involves the introduction of dUTP (the nucleotide deoxyuradine) in lieu of TTP (the nucleotide thymine), along with the enzyme uracil N-glycosylase (UNG, also known as uracil DNA glycosylase, UDG), into the qPCR master mix. As a result of this substitution, all amplified product will contain deoxyuradine instead of thymidine.

5.5.2 False negative prevention

False negatives can result from inhibitors from the environment, poor experimental design of the qPCR, poor primer design, the variability of the environmental processing, or contaminants like DNase and RNase that destroy the nucleic acid. There are several procedures that are reported to reduce inhibition, including the incorporation of BSA (bovine serum albumin), DMSO (dimethyl sulfoxide), betaine (N,N,N-trimethylglycine=[carboxymethyl]trimethylammonium), formamide, non-ionic detergents (e.g. Triton X-100, Tween 20 or Nonidet P-40), TMAC (tetramethylammonium chloride), 7-deaza-2'-deoxyguanosine (dC7GTP), or T4 gene 32 protein, to the PCR. However, there are no procedures that work with all applications.

6 DATA RECORDING, RECORD KEEPING, AND DATA EVALUATION

The Laboratory of Medical Zoology keeps all generated data on the server, so a third party can assess the results.

6.1 Data recording and record keeping

The Laboratory of Medical Zoology should retain copies of all data records for a minimum of five years after the project is completed. Electronically maintained data

should be backed up on a regular basis and stored in a separate location from the original data. All data recording should be checked by the laboratory and QA/QC managers for correctness and completeness, and each entry should be checked for accuracy of transcription.

6.1.1 Equipment

All equipment calibration and maintenance should be documented in appropriate logbooks. A schedule of equipment calibration and maintenance should be established with the procurement of each new piece of equipment, and a copy of the schedule should be posted in the laboratory near the equipment.

6.1.2 Reagents, kits, primer sets, and enzymes

Logbooks should be maintained for all reagents, kits, primer sets, and enzymes, and should document all pertinent information needed to identify possible sources of contamination, including the following:

- Product name or name of the primer
- Manufacturer
- Product number
- Formulations (reagents)
- Sequence (primers and probes)
- Receipt or preparation analysts' name and initials
- Storage location and location of components
- Location of the commercially available kit specifications and procedure
- Concentration of original primer stocks and working solutions
- Storage buffer
- Number of units (enzyme)
- Associated buffer (enzyme)

The reagent logbooks should be centrally located and readily accessible. See Appendix A for the policy of expired arthropod testing reagents.

6.1.3 Sample processing and analysis

The Laboratory of Medical Zoology documents all steps involved in the analysis of biological samples, including the handling, processing, and examination of the samples. Information (date, investigator, source, location) that identifies the samples are recorded, and the sample condition upon receipt. Detailed information should be recorded concerning the analysis of the sample; including the date and time each analytical step was performed, as well as the researcher. For all procedures, data should be presented in a clear and concise way to ensure that results are interpreted properly. All laboratory QA/QC procedures for the method should be followed. Before being reported to the end user, all data should be reviewed and eventually analyzed to assure that appropriate controls were used. Any information that impacts the quality or validity of the data should be noted on the data report.

6.1.4 Quantitative real-time PCR

Raw data, including fluorescence growth curves, the fluorescence threshold values, should be maintained. It should be saved and backed up by paper copy or in an electronic form other than the real-time instrument.

Appendix A



Policy for Expired Arthropod Testing Reagents August, 2017

1. Expired reagents are defined as those in which the manufacturer or certification can no longer guarantee quality.
2. No testing is permitted if reagents have exceeded their expiration date.
3. Label the arthropod testing reagents (e.g. master mix, extraction reagents) with the date it was opened or prepared and initial.
4. Monitor conditions where reagents are stored (e.g., temperature).
5. Discard expired reagents through UMass-Amherst EHS Program.